

Product Datasheet

Neuropeptide Y Antibody

NB100-1624

Unit Size: 0.05 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-1624

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-1624



NB100-1624

Neuropeptide Y Antibody

Product Information	
Unit Size	0.05 ml
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein G purified
Buffer	Whole serum


Product Description	
Description	Novus Biologicals Guinea Pig Neuropeptide Y Antibody (NB100-1624) is a polyclonal antibody validated for use in IHC and ICC/IF. Anti-Neuropeptide Y Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Guinea Pig
Gene ID	4852
Gene Symbol	NPY
Species	Human, Mouse, Rat
Reactivity Notes	Reacts with Human, Mouse and Rat.
Specificity/Sensitivity	Neuropeptide Y
Immunogen	Synthetic peptide: SDLLMRESTENAPRTR, corresponding to amino acids 76-91 of Rat Neuropeptide Y.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Immunohistochemistry 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500




Images

Immunohistochemistry: Neuropeptide Y Antibody [NB100-1624] - Image: NPY staining of striatal rat interneurons (arrows).



Immunohistochemistry: Neuropeptide Y Antibody [NB100-1624] - Hypothalamus



Publications

Varga A, Mészáros Z, Sivadón M et al. Spinal Excitatory Dynorphinergic Interneurons Contribute to Burn Injury-Induced Nociception Mediated by Phosphorylated Histone 3 at Serine 10 in Rodents International Journal of Molecular Sciences 2021-02-25 [PMID: 33669046] (Immunocytochemistry/ Immunofluorescence)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-1624

NB120-7134	Goat anti-Guinea Pig IgG (H+L) Secondary Antibody (Pre-adsorbed)
NBP1-97036-5mg	Guinea Pig IgG Isotype Control
H00004852-P02-10ug	Recombinant Human Neuropeptide Y GST (N-Term) Protein
DLP00	Leptin/OB [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-1624

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

