

# Product Datasheet

## AMPK alpha 2 Antibody - BSA Free NB100-238

Unit Size: 100 ul

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 7

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-238](http://www.novusbio.com/NB100-238)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-238](http://www.novusbio.com/reviews/destination/NB100-238)



**NB100-238**

AMPK alpha 2 Antibody - BSA Free


Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)
Target Molecular Weight	62 kDa

Product Description	
Description	Novus Biologicals Rabbit AMPK alpha 2 Antibody - BSA Free (NB100-238) is a polyclonal antibody validated for use in WB and IP. Anti-AMPK alpha 2 Antibody: Cited in 7 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5563
Gene Symbol	PRKAA2
Species	Human, Mouse
Reactivity Notes	Orangutan (100%).
Immunogen	The immunogen recognized by this antibody maps to a region between residues 350 and 400 of human 5'-AMP-activated Protein Kinase, catalytic Alpha-2 chain using the numbering given in entry AAH69823.1 (GeneID 5563).

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:5000, Immunoprecipitation 2-8 ug/mg lysate

**Images**

Western Blot: AMPK alpha 2 Antibody [NB100-238] - Detection of Human and Mouse AMPK alpha 2 by Western Blot. Samples: Whole cell lysate (50 ug) from HeLa, 293T, Jurkat, mouse Renca, and mouse TCMK-1 cells. Antibodies: Affinity purified rabbit anti-AMPK alpha 2 antibody NB100-238 used for WB at 0.4 ug/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.



## Publications

Masuda M, Miyazaki-Anzai S, Keenan AL et al. Saturated phosphatidic acids mediate saturated fatty acid-induced vascular calcification and lipotoxicity. *J Clin Invest* 2015-10-26 [PMID: 26517697]

Oral mucosa harbors a high frequency of endothelial cells - a novel postnatal cell source for angiogenic regeneration *Stem Cells Dev.*, 2016-12-22;0(0):. 2016-12-22 [PMID: 27832737]

Roan, J N, Cheng, H N Et al. Exendin-4, a glucagon-like peptide-1 analogue, accelerates diabetic wound healing. *J Surg Res* 2017-02-01 [PMID: 27993221]

Derkach K, Zakharova I, Zorina I et al. The evidence of metabolic-improving effect of metformin in Ay/a mice with genetically-induced melanocortin obesity and the contribution of hypothalamic mechanisms to this effect *PLoS ONE* 2019-03-14 [PMID: 30870482] (WB, Mouse)

Jaitovich A, Angulo M, Lecuona E et al. High CO<sub>2</sub> levels cause skeletal muscle atrophy via AMPK, FoxO3a and muscle-specific ring finger protein1 (MuRF1). *J Biol Chem.* 2015-02-17 [PMID: 25691571] (WB, Mouse)

### Details:

AMPK alpha 2 antibody used for WB application on total cell lysates as well as cytoplasmic fractions of C2C12 mouse myoblasts generated myotubes that were transfected with AMPK alpha 1 as well as AMPK alpha 2 siRNA followed by exposure to high CO<sub>2</sub> (FIGURE 8C, 9B).

Rubin LJ, Magliola L, Feng X et al. Metabolic activation of AMP kinase in vascular smooth muscle. *J Appl Physiol* 2005-01-01 [PMID: 15377643]

Gusarova GA, Dada LA, Kelly AM et al.  $\alpha$ 1-AMP-Activated Protein Kinase Regulates Hypoxia-Induced Na,K-ATPase Endocytosis via Direct Phosphorylation of Protein Kinase C $\zeta$ . *Mol Cell Biol* 2009;29(13):3455-3464. [PMID: 19380482]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-238**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-238](http://www.novusbio.com/reviews/submit/NB100-238)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

