

Product Datasheet

PUM1 Antibody - BSA Free

NB100-259

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 7

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-259

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-259



NB100-259

PUM1 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)
Target Molecular Weight	126 kDa


Product Description	
Description	Novus Biologicals Goat PUM1 Antibody - BSA Free (NB100-259) is a polyclonal antibody validated for use in WB, ICC/IF and IP. Anti-PUM1 Antibody: Cited in 6 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Goat
Gene ID	9698
Gene Symbol	PUM1
Species	Human, Mouse
Immunogen	The immunogen recognized by this antibody maps to a region between residues 225 and 275 of human Pumilio 1 using the numbering given in entry NP_055491.1 (GeneID 9698).

Product Application Details	
Applications	Western Blot, Gel Super Shift Assays, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000 - 1:10000, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation 2-10 ug/mg lysate, Gel Super Shift Assays 1:1-1:100
Application Notes	Gel Super Shift Assays and Immunoprecipitation were reported in scientific literature (PMID: 23071749). Use in ICC/IF was reported in scientific literature (PMID: 23071749).




Images


Western Blot: PUM1 Antibody [NB100-259] - Whole cell lysate (50 ug) from HEK293T, A-549, and K-562 cells prepared using NETN lysis buffer. Antibody: Affinity purified goat anti-Pumillio 1 antibody used for WB at 0.1 ug/ml. Secondary: HRPconjugated rabbit anti-goat IgG. Chemiluminescence with an exposure time of 10 seconds.



Western Blot: PUM1 Antibody [NB100-259] - Detection of mouse Pumillio 1 by western blot. Samples: Whole cell lysate (50 ug) from NIH 3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified goat anti-Pumillio 1 antibody NB100-259 used for WB at 0.1 ug/ml. Secondary: HRP-conjugated rabbit anti-goat IgG. Chemiluminescence with an exposure time of 10 seconds.



Immunoprecipitation: PUM1 Antibody [NB100-259] - Detection of human Pumillio 1 by western blot of immunoprecipitates. Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from HEK293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified goat anti-Pumillio 1 antibody NB100-259 (lot NB100-259-3) used for IP at 6 ug per reaction. Pumillio 1 was also immunoprecipitated by a previous lot of this antibody (lot NB100-259-2) and another rabbit anti-Pumillio 1 antibody. For blotting immunoprecipitated Pumillio 1, NB100-259 was used at 0.1 ug/ml. Chemiluminescence with an exposure time of 3 seconds.



Publications

Takei N, Takada Y, Kawamura S et al. Changes in subcellular structures and states of pumilio 1 regulate the translation of target Mad2 and cyclin B1 mRNAs *Journal of Cell Science* 2020-05-19 [PMID: 33148609]

Uyhazi KE, Yang Y, Liu N et al. Pumilio Proteins Exert Distinct Biological Functions and Multiple Modes of PostTranscriptional Regulation in Embryonic Stem Cell Pluripotency and Early Embryogenesis *bioRxiv* 2019-08-30

Ma L, Dong F, Zaid M et al. ABCA1 protein enhances Toll-like receptor 4 (TLR4)-stimulated interleukin-10 (IL-10) secretion through protein kinase A (PKA) activation. *J Biol Chem* 2012-11-23 [PMID: 23055522]

Vessey JP, Schoderboeck L, Gingl E et al. Mammalian Pumilio 2 regulates dendrite morphogenesis and synaptic function. *Proc Natl Acad Sci U S A* 2010-02-01 [PMID: 20133610]

Morris AR, Mukherjee N, Keene JD et al. Ribonomic analysis of human Pum1 reveals cis-trans conservation across species despite evolution of diverse mRNA target sets. *Mol Cell Biol* 2008-06-01 [PMID: 18411299]

Chiu YL, Witkowska HE, Hall SC et al. High-molecular-mass APOBEC3G complexes restrict Alu retrotransposition. *Proc Natl Acad Sci U S A* 2006-10-01 [PMID: 17030807]

Zhu X, Wang Y, Pi W et al. NF-Y recruits both transcription activator and repressor to modulate tissue- and developmental stage-specific expression of human gamma-globin gene. *PLoS One*. 2012-10-10 [PMID: 23071749] (GS, IP, Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-259

NB800-PC1	HeLa Whole Cell Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-259

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



