

Product Datasheet

CKII alpha prime polypeptide Antibody - BSA Free NB100-379

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 7

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-379

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-379



NB100-379**CKII alpha prime polypeptide Antibody - BSA Free**

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)
Target Molecular Weight	41 kDa


Product Description	
Description	Novus Biologicals Rabbit CKII alpha prime polypeptide Antibody - BSA Free (NB100-379) is a polyclonal antibody validated for use in WB and IP. Anti-CKII alpha prime polypeptide Antibody: Cited in 7 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1459
Gene Symbol	CSNK2A2
Species	Human, Mouse
Immunogen	The immunogen recognized by this antibody maps to a region between residue 350 and the C-terminus (residue 350) of human Casein Kinase 2, alpha 2 using the numbering given in entry NP_001887.1 (GeneID 1459).

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:200- 1:10000, Immunoprecipitation 2- 10 ug/mg lysate
Application Notes	By Western Blot, this antibody recognizes a band at ~40 kDa, representing CKII alpha'.




Images

Western Blot: CKII alpha prime polypeptide Antibody [NB100-379] - Detection of Human and Mouse CKII alpha' by Western Blot. Samples: Whole cell lysate (50 ug) from HeLa, 293T, Jurkat, mouse TCMK-1, and mouse NIH3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-CKII alpha' antibody NB100-379 used for WB at 0.1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Immunoprecipitation: CKII alpha prime polypeptide Antibody [NB100-379] - Detection of human CKII alpha' by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-CKII alpha' antibody NB100-379 used for IP at 6 ug per reaction. CKII alpha' was also immunoprecipitated by rabbit anti-CKII alpha' antibody NB100-378. For blotting immunoprecipitated CKII alpha', NB100-379 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 10 seconds.



Publications

Warsinger-Pepe N, Chang C, Desroberts CR et al. Polycomb response elements reduce leaky expression of Cas9 under temperature-inducible Hsp70Bb promoter in *Drosophila melanogaster* G3 (Bethesda) 2023-04-13 [PMID: 36705519]

Mansky RH, Greguske EA, Yu D et al. Tumor suppressor p53 regulates heat shock factor 1 protein degradation in Huntington's disease *Cell reports* 2023-03-02 [PMID: 36867535] (WB, Mouse)

Yu D, Zarate N, White A et al. CK2 alpha prime and alpha-synuclein pathogenic functional interaction mediates synaptic dysregulation in huntington's disease *Acta neuropathologica communications* 2022-06-03 [PMID: 35659303] (WB, Mouse)

Deshiere A, Duchemin-Pelletier E, Spreux E et al. Unbalanced expression of CK2 kinase subunits is sufficient to drive epithelial-to-mesenchymal transition by Snail1 induction. *Oncogene* 2012-05-01 [PMID: 22562247]

Hamdorf M, Berger A, Schule S et al. PKC Δ -induced PU.1 phosphorylation promotes hematopoietic stem cell differentiation to dendritic cells. *Stem Cells* 2011-02-01 [PMID: 21732487]

Gao Y, Wang HY. Casein kinase 2 is activated and essential for Wnt/beta-catenin signaling. *J Biol Chem* 2006-07-01 [PMID: 16672224]

Xu J, Reddy BJ, AnP et al. Casein kinase 2 reverses tail-independent inactivation of kinesin-1 *Nat Commun* 2012-03-27 [PMID: 22453827] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-379

NBL1-09537	CKII alpha prime polypeptide Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-379

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

