

Product Datasheet

PINK1 Antibody - BSA Free

NB100-493

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 12

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-493

Updated 1/19/2026 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-493



NB100-493

PINK1 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS
Target Molecular Weight	62.7 kDa

Product Description	
Description	Novus Biologicals Rabbit PINK1 Antibody - BSA Free (NB100-493) is a polyclonal antibody validated for use in WB. Anti-PINK1 Antibody: Cited in 12 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	65018
Gene Symbol	PINK1
Species	Human, Mouse, Rat
Reactivity Notes	Rat (PMID: 24411077) and mouse (PMID: 21760537) reactivity reported in scientific literature
Specificity/Sensitivity	PINK1 Antibody is expected recognize isoform 1 but will not recognize isoform 2.
Immunogen	PINK1 antibody was developed using an N-terminal region synthetic peptide made to the human PINK1 protein sequence (between residues 1-50). [UniProt# Q9BXM7]

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 2 ug/ml, Immunocytochemistry/ Immunofluorescence
Application Notes	This PINK1 antibody is useful for Western blot, where a band is seen ~63 kDa. The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors. Unprocessed PINK1 is 63 kDa which undergoes proteolytic processing to generate 55 kDa and 42 kDa cleaved forms. This antibody is made to a region that lies within the mitochondrial targeting sequence and is cleaved off to generate a mature protein.

Images

Western Blot: PINK1 Antibody [NB100-493] - Detection of murine PINK1 using NB100-493.


Lane 1: molecular weight marker.

Lane 2. MES cell Mitochondria (20 ug) with a band at the observed molecular weight of 63 kDa.

Lane 3. MES cytosol (20 ug).

Lane 4. MES nuclear (20 ug) as negative control.

Lane 5. Purified human cytochrome C (0.1 ug) as PINK1 negative control.



PINK1 was detected in immersion fixed HeLa human cervix adenocarcinoma cell line using Rabbit anti-PINK1 Antigen Affinity Purified Polyclonal Antibody (Catalog # NB100-493) at 1.0 $\mu\text{g}/\text{mL}$ overnight at 4C. Cells were stained using DyLight 488-conjugated Anti-Rabbit IgG (H+L) Cross-Absorbed Secondary Antibody (green), and counterstained with DAPI (blue). Cells were imaged using a 100X objective and digitally deconvolved.



Copyright © 2026 Novus Biologicals, a Bio-Techne Brand
Magnification: 100x

Publications

Hu WF, Lee CH, Pang CY et Al. Mitigating mitochondrial dysfunction and neuroinflammation by hematoma aspiration in a new surgical model for severe intracerebral hemorrhage *Exp Neurol* 2024-12-07 [PMID: 39653108]

Kaur, B;Miglioranza Scavuzzi, B;Yang, M;Yao, J;Jia, L;Abcouwer, SF;Zacks, DN; ER Stress and Mitochondrial Perturbations Regulate Cell Death in Retinal Detachment: Exploring the Role of HIF1? *Investigative ophthalmology & visual science* 2024-09-03 [PMID: 39325470]

Casey R Appell, Nigel C Jiwan, Rui Wang, Chwan-Li Shen, Hui-Ying Luk Ginger Supplementation Attenuated Mitochondrial Fusion and Improved Skeletal Muscle Size in Type 2 Diabetic Rats. In vivo (Athens, Greece) 2023-12-28 [PMID: 38148056]

Maheshwari C, Vidoni C, Titone R et al. Isolation, Characterization, and Autophagy Function of BECN1-Splicing Isoforms in Cancer Cells *Biomolecules* 2022-08-02 [PMID: 36008963] (WB, Human)

Chiang S, Braidy N, Maleki S et al. Mechanisms of impaired mitochondrial homeostasis and NAD⁺ metabolism in a model of mitochondrial heart disease exhibiting redox active iron accumulation *Redox Biol* 2021-08-20 [PMID: 34416478]

Ramasubramanian B, Griffith C, Hanson M et al. Protective Effects of Chaya against Mitochondrial and Synaptic Toxicities in the Type 2 Diabetes Mouse Model TallyHO Cells 2022-02-21 [PMID: 35203393] (WB, IF/IHC, Mouse)

Kim Sj, Cheresh P, Jablonski Rp et Al. Mitochondrial 8-Oxoguanine DNA Glycosylase Mitigates Alveolar Epithelial Cell PINK1 Deficiency, Mitochondrial DNA Damage, Apoptosis and Lung Fibrosis *Am. J. Physiol. Lung Cell Mol. Physiol.* 2020-03-25 [PMID: 32209025] (WB, Mouse)

Hwang, SH;Kim, MC;Ji, S;Yang, Y;Jeong, Y;Kim, Y; Glucose starvation induces resistance to metformin through the elevation of mitochondrial multidrug resistance protein 1 *Cancer Sci.* 2019-04-01 [PMID: 30689265] (WB, Human)

Amadoro G, Corsetti V, Florenzano F et al. AD-linked, toxic NH₂ human tau affects the quality control of mitochondria in neurons. *Neurobiol. Dis.* 2014-02-01 [PMID: 24411077] (WB, Rat)

Yasuda T, Hayakawa H, Nihira T et al. Parkin-mediated protection of dopaminergic neurons in a chronic MPTP-minipump mouse model of Parkinson disease *J Neuropathol Exp Neurol* 2011-08-01 [PMID: 21760537] (WB, Mouse)

Lutz AK, Exner N, Fett ME et al. Loss of Parkin or PINK1 Function Increases Drp1-dependent Mitochondrial Fragmentation. *J Biol Chem*;284(34):22938-22951. 2009-01-01 [PMID: 19546216]

Exner N, Treske B, Paquet D et al. Loss-of-function of human PINK1 results in mitochondrial pathology and can be rescued by parkin. *J Neurosci*;27(45):12413-8. 2007-11-07 [PMID: 17989306] (WB, Human)

Procedures

Western Blot protocol for PINK1 Antibody (NB100-493)

Western Blot Protocol

1. Perform SDS-PAGE on samples to be analyzed, loading 10-25 ug of total protein per lane.
2. Transfer proteins to PVDF membrane according to the instructions provided by the manufacturer of the membrane and transfer apparatus.
3. Stain the membrane with Ponceau S (or similar product) to assess transfer success, and mark molecular weight standards where appropriate.
4. Rinse the blot TBS -0.05% Tween 20 (TBST).
5. Block the membrane in 5% Non-fat milk in TBST (blocking buffer) for at least 1 hour.
6. Wash the membrane in TBST three times for 10 minutes each.
7. Dilute primary antibody in blocking buffer and incubate overnight at 4C with gentle rocking.
8. Wash the membrane in TBST three times for 10 minutes each.
9. Incubate the membrane in diluted HRP conjugated secondary antibody in blocking buffer (as per manufacturer's instructions) for 1 hour at room temperature.
10. Wash the blot in TBST three times for 10 minutes each (this step can be repeated as required to reduce background).
11. Apply the detection reagent of choice in accordance with the manufacturer's instructions.

Immunocytochemistry/ Immunofluorescence Protocol for PINK1 Antibody (NB100-493)

Immunocytochemistry Protocol

Culture cells to appropriate density in 35 mm culture dishes or 6-well plates.

1. Remove culture medium and wash the cells briefly in PBS. Add 4% paraformaldehyde to the dish and fix at room temperature for 10 minutes.
2. Remove the paraformaldehyde and wash the cells in PBS.
3. Permeabilize the cells with 0.1% Triton X100 or other suitable detergent for 2 min.
4. Remove the permeabilization buffer and wash three times for 5 minutes each in PBS. Be sure to not let the specimen dry out.
5. To block nonspecific antibody binding, incubate in 10% normal goat serum from 1 hour to overnight at room temperature.
6. Add primary antibody at appropriate dilution and incubate overnight at 4C.
7. Remove primary antibody and replace with PBS. Wash three times for 5 minutes each.
8. Add secondary antibody at appropriate dilution. Incubate for 1 hour at room temperature.
9. Remove secondary antibody and replace with PBS. Wash three times for 5 minutes each.
10. Counter stain DNA with DAPI if required.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-493

NB100-493PEP	PINK1 Antibody Blocking Peptide
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-493

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

