

Product Datasheet

ARID1A Antibody - BSA Free NB100-55334

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 5

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-55334

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-55334



NB100-55334

ARID1A Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)


Product Description	
Description	Novus Biologicals Rabbit ARID1A Antibody - BSA Free (NB100-55334) is a polyclonal antibody validated for use in WB and ChIP. Anti-ARID1A Antibody: Cited in 5 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8289
Gene Symbol	ARID1A
Species	Human, Mouse, Zebrafish
Reactivity Notes	Mouse, Zebrafish reactivity reported from verified customer reviews.
Immunogen	The immunogen recognized by this antibody maps to a region between residue 750 and 800 of human AT rich interactive domain 1A (BRG1-associated factor 250kD) using the numbering given in entry NP_006006.3 (GeneID 8289).

Product Application Details	
Applications	Western Blot, Chromatin Immunoprecipitation (ChIP), Immunoprecipitation (Negative)
Recommended Dilutions	Western Blot 1:2000-1:10000, Immunoprecipitation (Negative), Chromatin Immunoprecipitation (ChIP)
Application Notes	Use in chromatin immunoprecipitation reported in scientific literature (PMID: 30337622). Immunoprecipitation Not recommended. ARID1A antibody validated for WB from verified customer reviews.




Images


Western Blot: ARID1A Antibody [NB100-55334] - Whole cell lysate (5, 15 and 50 ug) from HeLa cells. Antibody used for WB at 0.1 ug/ml.




Chromatin Immunoprecipitation: ARID1A Antibody [NB100-55334] - In vivo Smarce1-binding sites were detected by chromatin immunoprecipitation using a SMARCE1 antibody. Cropped gel images representing typical results from ChIP assays in: the regulatory module located from the transcription initiation site +1 to -1415, which includes the basic promoter (*gata5-Bp*), the module located from -1469 to -1892 (*gata5-A*), the module located from -7823 to -8213 (*gata5-B*) and the module located from -16805 to -17147 (*gata5-C*). IgG, Smarce1, Smarca4, Arid1a and Pbrm indicate the PCR products from genomic DNA precipitated using the corresponding antibody. WT and Mut indicate the source of genomic DNA. The input is the PCR product from the total input genomic DNA. Full-length gels are presented in Supplementary Figure 3. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/s41598-018-33746-8>), licensed under a CC-BY license.



Western Blot: ARID1A Antibody [NB100-55334] - Mouse embryonic fibroblasts nuclear extract at 1:1000 dilution. (Nuclear extract instead of whole cell lysate provides much better signal.) Image provided via verified customer review.



Western Blot: ARID1A Antibody [NB100-55334] - Analysis of ARID1A in zebrafish embryos using anti-ARID1A antibody. Image from verified customer review.



Publications

Yakovlev VA, Sullivan SA, Fields EC, Temkin SM. PARP inhibitors in the treatment of ARID1A mutant ovarian clear cell cancer: PI3K/Akt1-dependent mechanism of synthetic lethality *Frontiers in Oncology* 2023-02-22 [PMID: 36910637]

Shi H, Tao T, Abraham B et al. ARID1A loss in neuroblastoma promotes the adrenergic-to-mesenchymal transition by regulating enhancer-mediated gene expression *Sci. Adv.* 2020-07-01 [PMID: 32832616] (WB, Zebrafish)

Bullard WL, Kara M, Gay LA et al. Identification of murine gammaherpesvirus 68 miRNA-mRNA hybrids reveals miRNA target conservation among gammaherpesviruses including host translation and protein modification machinery *PLoS Pathog.* 2019-08-01 [PMID: 31393953] (WB, Mouse)

Castillo-Robles J, Ramirez L, Spaink HP, Lomeli H. smarce1 mutants have a defective endocardium and an increased expression of cardiac transcription factors in zebrafish. *Sci Rep.* 2018-10-18 [PMID: 30337622] (Chemotaxis, Zebrafish)

Wurster AL, Precht P, Pazin MJ et al. NF- κ B and BRG1 bind a distal regulatory element in the IL-3/GM-CSF locus. *Mol Immunol* 2011-09-01 [PMID: 21831442]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-55334

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-55334

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

