

Product Datasheet

Caspase-9 Antibody - BSA Free NB100-56118

Unit Size: 0.05 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 28

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-56118

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-56118



NB100-56118

Caspase-9 Antibody - BSA Free


Product Information	
Unit Size	0.05 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Unpurified
Buffer	Whole antisera

Product Description	
Description	Novus Biologicals Rabbit Caspase-9 Antibody - BSA Free (NB100-56118) is a polyclonal antibody validated for use in IHC, WB, Flow, ICC/IF and IP. Anti-Caspase-9 Antibody: Cited in 28 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	842
Gene Symbol	CASP9
Species	Human, Mouse, Rat, Canine, Gerbil
Immunogen	Human Caspase 9 recombinant protein catalytic subunits (NP_001220).


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot reported in scientific literature (PMID 24184144), Immunohistochemistry, Immunocytochemistry/ Immunofluorescence reported in scientific literature (PMID 22223358), Immunoprecipitation 1:50-1:200, Immunohistochemistry-Paraffin 1:1000-1:5000, Immunohistochemistry-Frozen Assay Dependent, Flow (Intracellular)
Application Notes	Caspase-9 immunostaining may appear diffuse or punctate. It may be localized in the cytosol, mitochondria, and nucleus. Caspase-9 staining in the nucleus is considered to be an indication of active Caspase-9. Both pro and active/cleaved Caspase-9 staining may also be seen in the mitochondria. Caspase antibodies are classical tools for detecting inactive (pro) and active (cleaved) forms of the enzymes. The presence of the large or small subunits in western blots is considered to be a marker of caspase activation: the proform of caspase-9 is detected at ~50 kDa, the large subunit at ~35 kDa, and the small subunit at ~15 kDa. Intermediate caspase-9 cleavage forms may also be seen at ~21 kDa.

Images


Immunohistochemistry: Caspase-9 Antibody [NB100-56118] - Mice with two different spontaneous mutant alleles of Sharpin show variations in onset of CPDM phenotype. However, the onset is much more rapid in Sharpincpdm-Dem mice as evidenced by antibody detection of cleaved CASPASE 9 (T) when compared to Sharpincpdm mutants (R) or WT mice (Q,S). . Image collected and cropped by CiteAb from the following publication ([//doi.org/10.1371/journal.pone.0085666](https://doi.org/10.1371/journal.pone.0085666)) licensed under a CC-BY license.




Immunohistochemistry: Caspase-9 Antibody [NB100-56118] - Autophagy-deficient mice showed increased activity of apoptosis and necroptosis. Increased the expression of Caspase-9 in 12-week-old Atg7deltapan mice. Caspase-9 quantitation and representative IF microphotographs of Atg7F/F (n=5) and Atg7deltapan (n=5) pancreatic tissue stained for DAPI (blue) and Caspase-9 (green) (1/1000, scale bar=50?um. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/doi/finder/10.1038/cddis.2017.313>), licensed under a CC-BY license.




Immunohistochemistry-Paraffin: Caspase-9 Antibody [NB100-56118] - Analysis of Caspase-9 in mouse transplanted lung tumor section using anti-Caspase-9 antibody. Image from verified customer review.




Immunohistochemistry: Caspase-9 Antibody [NB100-56118] - Exacerbated apoptosis and altered inflammation in double-deficient Atg7deltapan-Rip3-/- mice. (c) Depletion of Rip3 and pancreatic Atg7 (Atg7deltapan-Rip3d/d) increased pancreatic caspase-9 in 12-week-old mice. Caspase-9 quantitation and representative IF microphotographs of Atg7deltapan (n=5) and Atg7deltapan-Rip3d/d (n=5) pancreatic tissue stained for DAPI (blue) and Caspase-9 (red) (1/1000).. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/doi/finder/10.1038/cddis.2017.313>), licensed under a CC-BY license.




Immunohistochemistry-Paraffin: Caspase-9 Antibody [NB100-56118] - Analysis of Dog Brain sections stained for Active/Cleaved Caspase-9 expression at 1:2000. A and B. The pattern of Caspase-9 staining may vary between different types of neurons. Hematoxylin-eosin counterstain.




Immunohistochemistry-Paraffin: Caspase-9 Antibody [NB100-56118] - Analysis of Dog Brain sections stained for cleaved Caspase-9 expression using this antibody at 1:2000. A) Section from a dog brain 2 hr after reperfusion injury. B) Section from a dog brain sham control (brain surgery but no injury). Hematoxylin-eosin counterstain.




Immunohistochemistry-Paraffin: Caspase-9 Antibody [NB100-56118] - Analysis Human Brain sections stained for cleaved Caspase-9 expression using this antibody at 1:2000. A) Section from a patient 24 hr after head trauma. B) Control: section from a patient with no known neurological disease or head injury. Hematoxylin-eosin counterstain.



Immunocytochemistry/ Immunofluorescence: Caspase-9 Antibody [NB100-56118] - Exacerbated apoptosis & altered inflammation in double-deficient *Atg7 Δ pan-Rip3 $^{-/-}$* mice. (a) Depletion of Rip3 & pancreatic Atg7 (*Atg7 Δ pan-Rip3d/d*) increased pancreatic caspase-3 in 12-week-old mice. Caspase-3 quantitation & representative IF microphotographs of *Atg7 Δ pan* (n=5) & *Atg7 Δ pan-Rip3d/d* (n=5) pancreatic tissue stained for DAPI (blue) & Caspase-3 (red) (anti-active caspase-3 ab2302, 1/50, scale bar=50 μ m). (b) Depletion of Rip3 & pancreatic Atg7 (*Atg7 Δ pan-Rip3d/d*) increased pancreatic Bax in 12-week-old mice. Bax quantitation & representative IF microphotographs of *Atg7 Δ pan* (n=5) & *Atg7 Δ pan-Rip3d/d* (n=5) pancreatic tissue stained for DAPI (blue) & Caspase-3 (red) (anti-Bax (sc-526, 1/50). (c) Depletion of Rip3 & pancreatic Atg7 (*Atg7 Δ pan-Rip3d/d*) increased pancreatic caspase-9 in 12-week-old mice. Caspase-9 quantitation & representative IF microphotographs of *Atg7 Δ pan* (n=5) & *Atg7 Δ pan-Rip3d/d* (n=5) pancreatic tissue stained for DAPI (blue) & Caspase-9 (red) (anti-caspase-9 NB100-56118, 1/1000). (d) Reduced infiltration inflammation of macrophages (F4/80, NBP2-12506, 1/75) & early T-lymphocytes (MPO, ab9535, 1/50) in double-deficient *Atg7 Δ pan-Rip3d/d* (n=10) mice compared with pancreatic *Atg7 Δ pan* (n=7). Data are mean \pm S.E.M. for the numbers of animals as indicated in the graph, *P<0.05, **P<0.01, NS = no significance Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/28703808>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Nuclear translocation of TCA cycle enzymes prevents doxorubicin-mediated cellular damage. a–c Enzyme activity of IDH-2 (n = 4) (a), MDH-2 (n = 3) (b), and PDHC (n = 3) (c) in nuclei isolated from cells transduced with tetracycline-inducible NLS-IDH-2 constructs. Activity is represented as the rate in nM/min or Δ OD/min, p value* <0.05, p value**** = 0.0001, and lysates were probed for PCNA to detect equal nuclear protein amounts used for the assay. d Reduced apoptotic cell death as measured by Annexin:V:FITC staining in cells expressing nuclear TCA cycle dehydrogenases prior to doxorubicin DRN treatment. N = 3 p value* <0.05 (e) Quantification of troponin I released from cardiomyocytes as a marker of cardiac injury, minimum n = 4 p value** <0.007, p value* <0.05. f Cell death as measured by Cell Titer GloTM in cells expressing nuclear dehydrogenases prior to DRN treatment. minimum n = 4 p value*** = 0.0001, p value**** <0.0001. NLS-EGFP-2A-NLS-mCherry was used as a control for the 2 A linker system. g DNA damage as assessed by γ -H2AX foci formation. Scale bar 20 μ m. h, i qRT-PCR showing the expression of anti-apoptotic genes BCL2 (h) and MCL1 (i). n = 3 P value ** <0.005, * <0.05. j–l Protein levels of anti-apoptotic proteins (j, k along with reduced caspase 9 cleavage (l) in cells expressing nuclear IDH-2 compared to cells treated with doxorubicin. Image collected and cropped by CiteAb from the following open publication (<https://pubmed.ncbi.nlm.nih.gov/37468519>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Zhou X, Xie L, Bergmann F et al. The bile acid receptor FXR attenuates acinar cell autophagy in chronic pancreatitis. *Cell Death Discov.* 2017-06-29 [PMID: 28660075] (Western Blot, Human)

Jaworski S, Strojny B, Sawosz E et al. Degradation of Mitochondria and Oxidative Stress as the Main Mechanism of Toxicity of Pristine Graphene on U87 Glioblastoma Cells and Tumors and HS-5 Cells *Int J Mol Sci* 2019-02-02 [PMID: 30717385] (Western Blot, Human)

Doke T, Huang S, Qiu C et al. Genome-wide association studies identify the role of caspase-9 in kidney disease *Science Advances* 2021-11-05 [PMID: 34739325] (Western Blot, Human)

Espinosa-Jeffrey A, Arrazola RA, Chu B et al. Trophic factors intervention regenerates the nestin-expressing cell population in a model of perinatal excitotoxicity: Implications for perinatal brain injury and prematurity *Integrative Molecular Medicine* 2016-06-25 [PMID: 35558521] (Western Blot, Human)

Yan Huo, Abudureheman Mijiti, Ruonan Cai, Zhaohua Gao, Maierpu Aini, Abudukadier Mijiti, Zhaoling Wang, Rui Qie Scutellarin alleviates type 2 diabetes (HFD/low dose STZ)-induced cardiac injury through modulation of oxidative stress, inflammation, apoptosis and fibrosis in mice. *Human & experimental toxicology* 2022-03-07 [PMID: 34610774]

Srivastava S, Gajwani P, Jousma J et al. Nuclear translocation of mitochondrial dehydrogenases as an adaptive cardioprotective mechanism *Nature communications* 2023-07-19 [PMID: 37468519] (WB, Human)

Details:

Dilutions: 1:1000

Feng Y, Yang L, Ma X et al. Photobiomodulation treatment inhibits neurotoxic astrocytic polarization and protects neurons in in vitro and in vivo stroke models *Neurochemistry international* 2022-12-17 [PMID: 36539162] (IHC-Fr, Mouse)

Li J, Camirand A, Zakikhani M et al. Parathyroid hormone related protein inhibition blocks triple negative breast cancer expansion in bone through epithelial to mesenchymal transition reversal *JBMR Plus* 2022-06-20 [PMID: 35720668]

Eskiler GG, Turna O, Ozkan AD et al. The response of the canine mammary simple carcinoma and carcinosarcoma cells to 5-aminolaevulinic acid-based photodynamic therapy: An in vitro study *Journal of photochemistry and photobiology. B, Biology* 2022-07-05 [PMID: 35810598] (WB, Canine)

Gao K, Niu J, Dang X Wnt-3a improves functional recovery through autophagy activation via inhibiting the mTOR signaling pathway after spinal cord injury *Neurosci. Lett.* 2020-08-17 [PMID: 32818590] (WB, Rat)

Xie L, Xia L, Klaiber U et al. Effects of neoadjuvant FOLFIRONOX and gemcitabine-based chemotherapy on cancer cell survival and death in patients with pancreatic ductal adenocarcinoma *Oncotarget* 2019-12-31 [PMID: 31921387] (ICC/IF, Human)

Zhou X, Xie L, Xia L et al. RIP3 attenuates the pancreatic damage induced by deletion of ATG7 *Cell Death and Disease* 2017-07-13 [PMID: 28703808] (ICC/IF, Mouse)

More publications at <http://www.novusbio.com/NB100-56118>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-56118

NBP3-11853	Jurkat Staurosporine Treated / Untreated Cell Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-56118

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

