

Product Datasheet

Histone H3 [p Ser28] Antibody (117C826) - BSA Free NB100-56595

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-56595

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-56595



NB100-56595

Histone H3 [p Ser28] Antibody (117C826) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	117C826
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	15 kDa


Product Description	
Description	Novus Biologicals Mouse Histone H3 [p Ser28] Antibody (117C826) - BSA Free (NB100-56595) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Histone H3 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	126961
Gene Symbol	H3C14
Species	Human, Fish
Reactivity Notes	Predicted to react with most Mammals and Invertebrates. Fish reactivity reported in scientific literature (PMID: 22279221).
Immunogen	This Histone H3 [p Ser28] antibody (117C826) was raised against a synthetic peptide containing phospho serine at position 28 CARKS*APATGGVKK of human Histone H3. The immunogen sequence is conserved across mammalian and invertebrates, as well as plants.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1-2 ug/ml, Immunocytochemistry/ Immunofluorescence 1:10-1:500. Use reported in scientific literature (PMID 22279221), Immunohistochemistry-Paraffin reported in scientific literature (PMID 35897708)
Application Notes	In anti-Fas mAb treated Jurkat, an 18 kDa band is observed using Histone H3 [p Ser28] antibody (117C826). Not recommended for Immunohistochemistry-Paraffin.



Images

Western Blot: Histone H3 [p Ser28] Antibody (117C826) [NB100-56595]
 - Analysis using Azide and BSA Free version of NB100-56595. Histone H3 phospho (Ser28) in (1) untreated, and (2) Fas antibody treated Jurkat lysate using Histone H3 phospho antibody at 2 ug/ml. Observed molecular weight is ~17 kDa



Publications

Montes-Mojarro IA, Hassas S, Staehle S et al. Multiparametric Classification of Non-Muscle Invasive Papillary Urothelial Neoplasms: Combining Morphological, Phenotypical, and Molecular Features for Improved Risk Stratification International journal of molecular sciences 2022-07-23 [PMID: 35897708] (IHC-P, Human)

Coffin AB, Mohr RA, Sisneros JA. Saccular-specific hair cell addition correlates with reproductive state-dependent changes in the auditory saccular sensitivity of a vocal fish. J Neurosci. 2012-01-25 [PMID: 22279221] (ICC/IF, Fish)

Details:

Histone H3 (Ser28)/PHH3 (IMG-371). IF/ICC: Fish epithelial cells, data described but not shown.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-56595

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-56595

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

