

Product Datasheet

PNUTS/PPP1R10 Antibody - BSA Free NB100-605

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 6

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-605

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-605



NB100-605

PNUTS/PPP1R10 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)


Product Description	
Description	Novus Biologicals Rabbit PNUTS/PPP1R10 Antibody - BSA Free (NB100-605) is a polyclonal antibody validated for use in WB and IP. Anti-PNUTS/PPP1R10 Antibody: Cited in 6 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5514
Gene Symbol	PPP1R10
Species	Human, Mouse
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID:31877321).
Immunogen	A synthetic peptide, mapping to a region region between residue 900 and the C-terminus (residue 940) of human Phosphatase Nuclear Targeting Subunit (a.k.a. Protein Phosphatase 1, regulatory subunit 10) using the numbering given in TrEMBL entry Q96QC0 (GeneID 5514).

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000-1:5000, Immunoprecipitation 1-10 ug/mg lysate
Application Notes	NB 100-605 may be used for Western blot. Suggested working dilutions: * Western Blot: 1:1,000 - 1:5,000 Immunoprecipitation: 1 to 10 ug/mg lysate Histo/Cytochemistry: not tested. *The investigator should determine the optimal working dilution for a specific application.




Images

Western Blot: PNUTS/PPP1R10 Antibody [NB100-605] - Detection of Human PNUTS by Western Blot. Samples: Whole cell lysate (50 ug) from HeLa and 293T cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-PNUTS antibody NB100-605 used for WB at 0.2 ug/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.



Western Blot: PNUTS/PPP1R10 Antibody [NB100-605] - Detection of Human PNUTS on HeLa whole cell lysates using NB100-605. IP of PNUTS was successfully performed with NB100-603, NB100-604, and another PNUTS polyclonal. IPed PNUTS was blotted using NB100-605



Publications

Travis B McMurphy Hepatic Expression of Adenovirus 36 E4ORF1 Improves Glycemic Control and Promotes Glucose Metabolism via AKT Activation Diabetes, 2016-11-30;0(0):. 2016-11-30 [PMID: 27903748]

F Nakazeki, M Nishiga, T Horie, H Nishi, Y Nakashima, O Baba, Y Kuwabara, T Nishino, T Nakao, Y Ide, S Koyama, M Kimura, S Tsuji, N Sowa, S Yoshida, SJ Conway, M Yanagita, T Kimura, K Ono Loss of periostin ameliorates adipose tissue inflammation and fibrosis in vivo Sci Rep, 2018-06-04;8(1):8553. 2018-06-04 [PMID: 29867212]

Kwon J, Yeh Y, Kawarasaki S et al. Mevalonate biosynthesis pathway regulates the development and survival of brown adipocytes iScience 2023-03-01 [PMID: 36895651]

Pathak, P;Chiang, JYL; Sterol 12 alpha-hydroxylase Aggravates Dyslipidemia by Activating the Ceramide/mTORC1/SREBP1C Pathway via FGF21 and FGF15 Gene Expr. 2019-03-19 [PMID: 30890204]

Jing X, Wu J, Dong C et al. COVID-19 instigates adipose browning and atrophy through VEGF in small mammals Nature metabolism 2022-12-01 [PMID: 36482111]

Panagiotakopoulou V, Botsakis K, Delis F et al. Anti-neuroinflammatory, protective effects of the synthetic microneurotrophin BNN-20 in the advanced dopaminergic neurodegeneration of "weaver" mice Neuropharmacology 2020-03-15 [PMID: 31877321] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-605

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-605

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

