

Product Datasheet

INTS11 Antibody - BSA Free

NB100-60638

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-60638

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-60638



NB100-60638

INTS11 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)
Target Molecular Weight	68 kDa


Product Description	
Description	Novus Biologicals Rabbit INTS11 Antibody - BSA Free (NB100-60638) is a polyclonal antibody validated for use in IHC, WB and IP. Anti-INTS11 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	54973
Gene Symbol	INTS11
Species	Human, Mouse
Immunogen	The immunogen recognized by this antibody maps to a region between residue 550 and 600 of human integrator complex subunit 11 using the numbering given in entry NP_060341.2 (GenelD 54973).

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000-1:10000, Immunohistochemistry 1:1000 - 1:5000, Immunoprecipitation 2-10 ug/mg lysate, Immunohistochemistry-Paraffin 1:1000 - 1:5000
Application Notes	Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections. INTS11 antibody validated for WB from a verified customer review.




Images


Western Blot: INTS11 Antibody [NB100-60638] - Whole cell lysate (50 ug) from HeLa, K-562, A-549, TCMK-1, and NIH 3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-INT11 antibody used for WB at 0.1 ug/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.



Immunohistochemistry-Paraffin: INTS11 Antibody [NB100-60638] - Section of human breast carcinoma. Antibody: Affinity purified rabbit anti-INT11 used at a dilution of 1:5,000 (0.2ug/ml). Detection: DAB




Western Blot: INTS11 Antibody [NB100-60638] - analysis of INTS11 in HeLa whole cell extracts (samples treated with siCtrl or siINTS11) using anti-INTS11 antibody. Image from verified customer review.



Whole cell extracts from HeLa cells treated with siCtrl or siINTS11 were resolved in SDS-PAGE.

Immunoprecipitation: INTS11 Antibody [NB100-60638] - Detection of human INT11 by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-INT11 antibody NB100-60638 (lot NB100-60638-3) used for IP at 6 ug per reaction. INT11 was also immunoprecipitated by a previous lot of this antibody (lot NB100-60638) and another rabbit anti-INT11 antibody. For blotting immunoprecipitated INT11, NB100-60638 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Publications

Cihlarova Z, Kubovciak J, Sobol M et al. BRAT1 links Integrator and defective RNA processing with neurodegeneration Nature communications 2022-08-26 [PMID: 36028512] (WB, Human)

Details:

Dilutions: 1:1000





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-60638

NBL1-09447	INTS11 Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-60638

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

