

Product Datasheet

Cystatin C Antibody (Cyst24cc) - BSA Free NB100-62127

Unit Size: 0.2 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-62127

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-62127




NB100-62127

Cystatin C Antibody (Cyst24cc) - BSA Free


Product Information	
Unit Size	0.2 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	Cyst24cc
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Product Description	
Description	Novus Biologicals Mouse Cystatin C Antibody (Cyst24cc) - BSA Free (NB100-62127) is a monoclonal antibody validated for use in ELISA. Anti-Cystatin C Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	1471
Gene Symbol	CST3
Species	Human
Immunogen	Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with Cystatin C purified from human urine.
Product Application Details	
Applications	ELISA, Immunoassay, Sandwich ELISA
Recommended Dilutions	ELISA, Immunoassay, Sandwich ELISA
Application Notes	Recommended pairs for sandwich immunoassay (capture-detection): Cyst24 - Cyst19, (NB100-62127/NBP1-05766) Cyst24 - Cyst28, (NB100-62127/NBP1-05765) Cyst23 - Cyst13, (NB100-62126/NBP1-05763)

Images

ELISA: Cystatin C Antibody (Cyst24) [NB100-62127] - Calibration curve
 One-step fluoroimmunoassay in streptavidin coated plates. Capture MAbs: Cyst24 and Cyst23 are biotinylated (200 ng/well). Detection MAbs: Cyst10, Cyst19, Cyst28 or Cyst13 are Eu3+- labelled (200 ng/ml). Incubation volume 100 μ l.



Sandwich ELISA: Cystatin C Antibody (Cyst24) [NB100-62127] - Titration curve of pooled normal serum in Cyst24. Cyst10- Eu3plus sandwich fluoroimmunoassay.



Publications

Thiriet PE, Medagoda D, Porro G, Guiducci C Rapid Multianalyte Microfluidic Homogeneous Immunoassay on Electrokinetically Driven Beads Biosensors 2020-12-21 [PMID: 33371213]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-62127

NBL1-09546	Cystatin C Overexpression Lysate
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-62127

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

