

Product Datasheet

LRP-1 Antibody (A2Mr alpha-2) - BSA Free NB100-64808

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-64808

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-64808



NB100-64808


LRP-1 Antibody (A2Mr alpha-2) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	A2Mr alpha-2
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Mouse LRP-1 Antibody (A2Mr alpha-2) - BSA Free (NB100-64808) is a monoclonal antibody validated for use in IHC, WB, ELISA, Flow, ICC/IF and IP. Anti-LRP-1 Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	4035
Gene Symbol	LRP1
Species	Human
Specificity/Sensitivity	Recognises an epitope in the extracellular region of the 500kD alpha chain.
Immunogen	Purified alpha2 macroglobulin receptor
Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot Non-reducing conditions, Flow Cytometry 1:50-1:100, ELISA Direct, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500
Application Notes	For Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Use non-reducing conditions for western blot.



Images

Immunohistochemistry: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Immunoperoxidase staining of a human tonsil cryosection with Mouse anti Human CD91 antibody, clone A2Mr alpha-2 followed by the Histar detection system. High power.




Immunohistochemistry: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Immunoperoxidase staining of a human tonsil cryosection with Mouse anti Human CD91 antibody, clone A2Mr alpha-2 followed by the Histar detection system. Low power.




Immunohistochemistry: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Immunoperoxidase staining of a human tonsil cryosection with Mouse anti Human CD91 antibody, clone A2Mr alpha-2 followed by the Histar detection system. Medium power.




Immunohistochemistry: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Immunoperoxidase staining of a human tonsil cryosection with Mouse anti Human CD91 antibody, clone A2Mr alpha-2 followed by the Histar detection system. Medium power.




Immunofluorescence: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Staining of a human tonsil cryosection with Mouse anti Human CD91 antibody, clone A2Mr alpha-2, red in A and Mouse anti Human CD21, clone LB21, green in B. C is the merged image with nuclei counterstained blue using DAPI. Low power.




Immunofluorescence: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Staining of a human tonsil cryosection with Mouse anti Human CD91 antibody, clone A2Mr alpha-2, red in A and Mouse anti Human CD21, clone LB21, green in B. C is the merged image with nuclei counterstained blue using DAPI. Medium power.




Immunofluorescence: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Staining of a human tonsil cryosection with Mouse anti Human CD91 antibody, clone A2Mr alpha-2, red in A and Mouse anti Human CD21, clone LB21, green in B. C is the merged image with nuclei counterstained blue using DAPI. High power.



Immunocytochemistry/ Immunofluorescence: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Reduced EHD1 recruitment to endosomes upon EHD4 knock-down. A-F, Representative micrographs & insets depicting EHD1-GFP recruitment to endosomes in untreated (A & inset in D), mock-treated (B & inset in E), & EHD4 knock-down (C & inset in F) cells. CRISPR/Cas9 gene-edited NIH3T3 cells expressing endogenous levels of EHD1 with GFP fused to the C-terminus (EHD1-GFP) were either untreated, mock-treated with transfection reagent, or transfected with EHD4 siRNA for 72 h. Cells were then incubated with anti-LRP1 antibody (30 min on ice, 30 min at 37°C), fixed, & imaged via confocal microscopy. G, Immunoblot showing reduced EHD4 (but not EHD1-GFP) expression in EHD1-GFP cells, with actin used as a loading control. The nitrocellulose filter paper was then stripped & immunoblotted with anti-GFP to show EHD1-GFP expression upon EHD4 loss. H, Graph depicting the mean count of EHD1-labeled endosomes in untreated, mock-treated & EHD4-depleted cells. Individual experiments were performed 3 times. Error bars denote standard deviation & p-values were determined by one-way ANOVA for individual experiments using a post-hoc Tukey HSD calculator to determine significance. A consensus p-value was then derived as described in the Materials & methods to assess significant differences between samples from the 3 experiments. Micrographs are representative orthogonal projects from three independent experiments, with 10 sets of z-stacks collected for each treatment per experiment. Bar, 10 μ m. n.s., not significant (consensus $p > 0.5$). Consensus p-values from Tukey HSD: * $p < 0.00001$. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32966336>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Immunocytochemistry/ Immunofluorescence: LRP-1 Antibody (A2Mr alpha-2) [NB100-64808] - Increased sorting endosome size & decreased EHD1 recruitment upon Rabenosyn-5 knock-down. A-L, Representative micrographs & insets depicting EEA1-labeled endosomes & EHD1-GFP in mock-treated & Rabenosyn-5 knock-down cells. CRISPR/Cas9 gene-edited NIH3T3 EHD1-GFP cells were either mock-treated with transfection reagent (A-F) or treated with Rabenosyn-5 siRNA (G-L) for 72 h. Cells were then incubated with anti-LRP1 antibody (30 min on ice, 30 min at 37°C), fixed & immunostained using anti-EEA1, & imaged by confocal microscopy. M, Immunoblot showing reduced Rabenosyn-5 expression in EHD1-GFP NIH3T3 cells. N, Graph depicting mean endosome size of mock-treated & Rabenosyn-5 knock-down cells. O, Graph depicting EHD1 recruitment to endosomes in mock-treated & Rabenosyn-5 knock-down cells. Error bars denote standard deviation & p-values were determined by independent two-tailed t-test, with significance derived from consensus p-values from the 3 experiments. Micrographs are representative orthogonal projections from three independent experiments, with 10 sets of z-stacks collected for each treatment per experiment. Bar, 10 μ m. ** $p < 0.00001$. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32966336>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Yashchenko A, Bland SJ, Song CJ et al. Cx3cr1 controls kidney resident macrophage heterogeneity *Frontiers in immunology* 2023-05-15 [PMID: 37256130] (FLOW, Mouse)

Jones T, Naslavsky N, Caplan S Eps15 Homology Domain Protein 4 (EHD4) is required for Eps15 Homology Domain Protein 1 (EHD1)-mediated endosomal recruitment and fission *PLoS ONE* 2020-09-23 [PMID: 32966336] (KD, ICC/IF, Mouse)

Dhawan K, Naslavsky N, Caplan S. Sorting nexin 17 (SNX17) links endosomal sorting to Eps15 Homology Domain Protein 1 (EHD1)-mediated fission machinery *J. Biol. Chem.* 2020-02-10 [PMID: 32041776]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-64808

NBL1-12668	LRP-1 Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-64808

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

