

Product Datasheet

PLK1 Antibody (13E8) NB100-74502

Unit Size: 100ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-74502

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-74502



NB100-74502


PLK1 Antibody (13E8)

Product Information	
Unit Size	100ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	13E8
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Purity	Protein G purified
Buffer	PBS with 1 mg/ml BSA
Product Description	
Description	Novus Biologicals Mouse PLK1 Antibody (13E8) (NB100-74502) is a monoclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-PLK1 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5347
Gene Symbol	PLK1
Species	Human, Mouse, Rat
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information.
Marker	Mitosis Marker
Specificity/Sensitivity	Detects PLK1 from human, mouse, and rat samples.
Immunogen	Recombinant protein containing residues 300-400 of human PLK1.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:20, Immunocytochemistry/Immunofluorescence 1:20, Immunoprecipitation 1:10 - 1:500, Immunohistochemistry-Paraffin 1:20




Images


Western Blot: PLK1 Antibody (13E8) [NB100-74502] - Proteins were transferred to a PVDF membrane and blocked with 5% Milk/TBST for at least 1 hour. Membranes were probed with PLK1 mouse monoclonal antibody at a dilution of 1:1000 overnight at 4C on a rocking platform. Membranes were washed in TBS-0.1%Tween 20 and probed with a goat anti-mouse-HRP secondary antibody at a dilution of 1:20,000 for at least one hour. Membranes were washed and chemiluminescent detection performed.




Immunocytochemistry/Immunofluorescence: PLK1 Antibody (13E8) [NB100-74502] - PLK1 staining (green), F-Actin staining with Phalloidin (red) and nuclei with DAPI (blue) in WiDr colon carcinoma cells. Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with or an antibody recognizing PLK1 at a dilution of 1:20 over night at 4C, washed with PBS and incubated with a DyLight-488 conjugated.




Immunohistochemistry-Paraffin: PLK1 Antibody (13E8) [NB100-74502] - Both normal and cancer biopsies of deparaffinized human Colon tissue.



Immunocytochemistry/Immunofluorescence: PLK1 Antibody (13E8) [NB100-74502] - PLK1 staining (green), F-Actin staining with Phalloidin (red) and nuclei with DAPI (blue) in HeLa cells. Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with or an antibody recognizing PLK1 at a dilution of 1:20 over night at 4C, washed with PBS and incubated with a DyLight-488 conjugated.



Immunocytochemistry/Immunofluorescence: PLK1 Antibody (13E8) [NB100-74502] - PLK1 staining (green), F-Actin staining with Phalloidin (red) and nuclei with DAPI (blue) in U251 glioma cells. Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with or an antibody recognizing PLK1 at a dilution of 1:20 over night at 4C, washed with PBS and incubated with a DyLight-488 conjugated.



Publications

Lee DS, Lee S, Kim C et al. The Casein Kinase 2 Inhibitor CX-4945 Promotes Cholangiocarcinoma Cell Death Through PLK1 Anticancer research 2022-07-01 [PMID: 35790285]

Chow M, Gu Y, He L Et al. Prognostic and Therapeutic Potential of The OIP5 Network in Papillary Renal Cell Carcinoma Research Square 2021-07-30 [PMID: 34503297]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-74502

NBL1-14515	PLK1 Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-74502

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

