

Product Datasheet

UTP18 Antibody - BSA Free

NB100-77295

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-77295

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-77295



NB100-77295

UTP18 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)


Product Description	
Description	Novus Biologicals Rabbit UTP18 Antibody - BSA Free (NB100-77295) is a polyclonal antibody validated for use in IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	51096
Gene Symbol	UTP18
Species	Human
Immunogen	The immunogen recognized by this antibody maps to a region between residue 325 and 375 of human small subunit (SSU) processome component UTP18 (WD repeat domain 50) using the numbering given in entry NP_057085.2 (GeneID 51096).

Product Application Details	
Applications	Immunoprecipitation
Recommended Dilutions	Immunoprecipitation 1:10-1:500
Application Notes	This antibody is not recommended for Western Blot analysis.




Images

Immunoprecipitation: UTP18 Antibody [NB100-77295] - Detection of Human UTP18/WDR50 by Western Blot of Immunoprecipitates. Samples: Whole cell lysate (1 mg for IP, 20% of IP loaded) from HeLa cells. Antibodies: Affinity purified rabbit anti-UTP18/WDR50 antibody NB100-77295 used for IP at 3 ug/mg lysate. UTP18/WDR50 was also immunoprecipitated by rabbit anti-UTP18/WDR50 antibody NB100-77294, which recognizes an upstream epitope. Detection: Chemiluminescence with an exposure time of 1 second.



Immunoprecipitation: UTP18 Antibody [NB100-77295] - Detection of human UTP18/WDR50 by western blot of immunoprecipitates. Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from HEK293T cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-UTP18/WDR50 antibody NB100-77295 used for IP at 3 ug per reaction. UTP18/WDR50 was also immunoprecipitated by rabbit anti-UTP18/WDR50 antibody NB100-77294. For blotting immunoprecipitated UTP18/WDR50, NB100-77294 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-77295

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-13511PEP	UTP18 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-77295

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

