

Product Datasheet

CD14 Antibody (M5E2) - BSA Free NB100-77758

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 14

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-77758

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-77758




NB100-77758

CD14 Antibody (M5E2) - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	M5E2
Preservative	0.02% Sodium Azide
Isotype	IgG2a Kappa
Purity	Protein G purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Mouse CD14 Antibody (M5E2) - BSA Free (NB100-77758) is a monoclonal antibody validated for use in IHC, WB, Flow, Dual RNAscope ISH-IHC, ICC/IF and Simple Western. Anti-CD14 Antibody: Cited in 14 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	929
Gene Symbol	CD14
Species	Human, Porcine, Bovine, Canine, Primate
Reactivity Notes	Expected to cross react with Chimpanzee, Cynomolgus, Rhesus, Pigtailed Macaque, Capuchin Monkey, Cotton-topped Tamarin, Common Marmoset, Squirrel Monkey.
Immunogen	Human peripheral blood cells.
Product Application Details	
Applications	Simple Western, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Block/Neutralize, CyTOF-ready, Dual RNAscope ISH-IHC
Recommended Dilutions	Simple Western 50 ug/ml, Flow Cytometry < 2.0 ug/10 ⁶ cells in 100 ul, Immunohistochemistry 1:10 - 1:500, Immunocytochemistry/ Immunofluorescence 1:10 - 1:500, Immunohistochemistry-Frozen 1:10 - 1:500, CyTOF-ready, Block/Neutralize reported in scientific literature (PMID 27228163), Dual RNAscope ISH-IHC
Application Notes	The M5E2 antibody inhibits monocyte activation and cytokine production induced by LPS. Additional reported applications (for the relevant formats) include: immunohistochemical staining of acetone-fixed frozen sections, blocking of LPS stimulation, and immunofluorescence microscopy. Clone M5E2 is not recommended for immunohistochemical staining of formalin-fixed paraffin-embedded sections. This antibody is CyTOF ready. See Simple Western Antibody Database for Simple Western validation: tested in THP-1; antibody dilution of 50 ug/ml; separated by size; detects a band at 54 kDa

Images


Immunocytochemistry/Immunofluorescence: CD14 Antibody (M5E2) [NB100-77758] - Representative immunostaining of CD14 (red) in adult human lung tissue sample with nuclei counterstained blue. This image was submitted by customer review. Image from the Alexa Fluor 647 version of this antibody.




Flow Cytometry: CD14 Antibody (M5E2) [NB100-77758] - A surface stain was performed on human peripheral blood monocytes with CD163 (EDHu-1) antibody NB110-40686APC and a matched isotype control NBP2-27287APC. Cells were incubated in an antibody dilution of 1 ug/mL for 20 minutes at room temperature. A co-stain was performed with NB100-77758AF488.




Immunocytochemistry/Immunofluorescence: CD14 Antibody (M5E2) [NB100-77758] - Human lung tissue was stained with FITC-conjugated CD14 (Green) for overnight. Image from the Alexa Fluor 488 version of this antibody.




Flow Cytometry: CD14 Antibody (M5E2) [NB100-77758] - Human peripheral blood monocytes stained with M5E2 APC




Flow Cytometry: CD14 Antibody (M5E2) [NB100-77758] - Analysis using the FITC conjugate of NB100-77758. Staining of Human peripheral blood monocytes with M5E2 FITC.




Flow Cytometry: CD14 Antibody (M5E2) [NB100-77758] - Analysis using the PE conjugate of NB100-77758. Staining of CD14 in human PBMC using anti-CD14 antibody conjugated with PE. Image from verified customer review.



Flow Cytometry: CD14 Antibody (M5E2) [NB100-77758] - Analysis using PE conjugate of NB100-77758. A cell surface stain was performed on hPBMCs with ILT5 antibody (MM0413-9S32) NBP2-11729 (top image) and a matched isotype control NBP2-27287 (bottom image). Cells were incubated in an antibody dilution of 1:200 for 20 minutes at room temperature. Both antibodies were conjugated to APC. A co-stain was performed using CD14 antibody (M5E2) NB100-77758PE.



Dual RNAscope ISH-IHC: CD14 Antibody (M5E2) [NB100-77758] - TNFA mRNA (red) and CD14 protein (green) were detected in formalin-fixed paraffin-embedded tissue sections of human malignant lymph node. ACD's Integrated Co-Detection Workflow was performed using ACD RNAscope Probe Hs-TNFA and CD14 antibody (M5E2) at 1:200 dilution. Tissue was stained on Leica Bond RX using RNAscope (TM) 2.5 LS Reagent Kit-RED, BOND Polymer Refine Detection (DAB) and Hematoxylin, BOND Polymer Refine Red Detection and Hematoxylin and RNAscope (TM) 2.5 LS Green Accessory Pack. Tissue was counterstained with 50% hematoxylin (blue).



Publications

Spindler LM, Feuerhake A, Ladel S et al. Nano-in-Micro-Particles Consisting of PLGA Nanoparticles Embedded in Chitosan Microparticles via Spray-Drying Enhances Their Uptake in the Olfactory Mucosa Frontiers in Pharmacology 2021-09-01 [PMID: 34539414]

Guido Trentadue, Gursah Kats-Ugurly, Tjasso Blokzijl, Gilles FH Diercks, Jan Willem Haveman, Klaas Nico Faber et al. Safe and Successful Treatment of Acute Cellular Rejection of an Intestine and Abdominal Wall Transplant With Vedolizumab Transplantation Direct 2020-02-01 [PMID: 32095513]

Kakiuchi-Kiyota, S, Obert, L A Et al. Expression of Hematopoietic Stem and Endothelial Cell Markers in Canine Hemangiosarcoma. Toxicol Pathol 2020-04-01 [PMID: 31918642] (IF/IHC, WB, Human)

Nitta S, Kusakari Y, Yamada Y et al. Conversion of mesenchymal stem cells into a canine hepatocyte-like cells by Foxa1 and Hnf4a Regenerative Therapy 2020-01-01 [PMID: 32123700] (FLOW, Canine)

Amin SB, Anderson KJ, Boudreau CE et al. Comparative Molecular Life History of Spontaneous Canine and Human Gliomas Cancer Cell 2020-02-10 [PMID: 32049048]

Bertolo A, Steffen F, Malonzo-Marty C et al Canine Mesenchymal Stem Cell Potential and the Importance of Dog Breed: Implication for Cell-Based Therapies. Cell Transplant. 2014-11-05 [PMID: 25375819] (FLOW, Canine)

Details:

Citation using the FITC version of this antibody.

Bertolo A, Pavlicek D, Gemperli A et al. Increased motility of mesenchymal stem cells is correlated with inhibition of stimulated peripheral blood mononuclear cells in vitro. J Stem Cells Regen Med 2017-12-18 [PMID: 29391751]

Details:

This citation used the FITC version of this antibody.

Walksman R, Lipinski MJ, Acampado E et al. Comparison of Acute Thrombogenicity for Metallic and Polymeric Bioabsorbable Scaffolds Magmaris Versus Absorb in a Porcine Arteriovenous Shunt Model. Circ Cardiovasc Interv. 2017-08-10 [PMID: 28801538] (Porcine)

Bertolo A, Capossela S, Frankl G et al. Oxidative status predicts quality in human mesenchymal stem cells Stem Cell Res Ther Jan 6 2017 12:00AM [PMID: 28061861] (FLOW, Human). Stem Cell Res Ther. 2017-01-06 [PMID: 28061861] (FLOW, Human)

Details:

This citation used the FITC version of this antibody.

Xun C Brown SA Jennings CD et al. Acute graft-versus-host-like disease induced by transplantation of human activated natural killer cells into SCID mice. Transplantation. [PMID: 8356598]

Rhoades JL Cibull ML Thompson JS et al. Role of natural killer cells in the pathogenesis of human acute graft-versus-host disease. Transplantation. [PMID: 8333033]

Waid TH Lucas BA Thompson JS et al. Treatment of acute cellular rejection with T10B9.1A-31 or OKT3 in renal allograft recipients. Transplantation. [PMID: 1531095]

More publications at <http://www.novusbio.com/NB100-77758>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-77758

NBL1-08905	CD14 Overexpression Lysate
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-77758

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

