

Product Datasheet

BANP Antibody - BSA Free NB100-79784

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-79784

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-79784



NB100-79784

BANP Antibody - BSA Free


Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)

Product Description	
Description	Novus Biologicals Rabbit BANP Antibody - BSA Free (NB100-79784) is a polyclonal antibody validated for use in IHC, WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	54971
Gene Symbol	BANP
Species	Human, Mouse, Rat
Immunogen	The immunogen recognized by this antibody maps to a region between residue 425 and the C-terminus (residue 469) of human BTG3 associated nuclear protein using the numbering given in TrEMBL entry Q96GJ7 (GeneID 54971).


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:10000-1:25000, Immunohistochemistry 1:10-1:500, Immunoprecipitation 1:500-1:2500, Immunohistochemistry-Paraffin 1:10-1:500
Application Notes	Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissues.

Images


Western Blot: BANP Antibody [NB100-79784] - Detection of human and mouse SMAR1 by western blot. Samples: Whole cell lysate (50 ug) from HeLa, HEK293T, and mouse NIH 3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-SMAR1 antibody NB100-79784 used for WB at 0.1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.




Immunohistochemistry: BANP Antibody [NB100-79784] - Detection of human SMAR1 by immunohistochemistry. Samples: FFPE sections of human lung carcinoma. Antibodies: Affinity purified rabbit anti-SMAR1 Cat. No. NB100-79784 Lot 2 used at a dilution of 1:5,000 (0.2ug/ml) (left) and Cat. No. NB100-79784 Lot3 used at a dilution of 1:5,000 (0.2ug/ml) (right). Detection: DAB



Western Blot: BANP Antibody [NB100-79784] - Nuclear extract (NE) from HeLa cells. Antibody used at 0.1 ug/ml.




Western Blot: BANP Antibody [NB100-79784] - Detection of Human, Mouse and Rat SMAR1/BANP by Western Blot. Samples: Whole cell lysate (50 ug) from mouse 4T1, mouse CT26.WT, rat C6, HeLa, and 293T cells. Antibodies: Affinity purified rabbit anti-SMAR1/BANP antibody NB100-79784 used for WB at 0.1 ug/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.




Immunohistochemistry-Paraffin: BANP Antibody [NB100-79784] - FFPE section of human breast tumor. Affinity purified rabbit anti-SMAR1/BANP used at a dilution of 1:250. Detection: DAB staining using Immunohistochemistry Accessory Kit




Immunohistochemistry: BANP Antibody [NB100-79784] - Sample: FFPE sections of human lung carcinoma (left) and mouse renal cell carcinoma (right). Antibody: Affinity purified rabbit anti-fs18 SMAR1 fs20 used at a dilution of 1:5,000 (0.2ug/ml). Detection: DAB



Immunohistochemistry: BANP Antibody [NB100-79784] - Detection of mouse SMAR1 by immunohistochemistry. Samples: FFPE sections of mouse renal cell carcinoma. Antibodies: Affinity purified rabbit anti-SMAR1 Cat. No. NB100-79784 Lot 2 used at a dilution of 1:5,000 (0.2ug/ml) (left) and Cat. No. NB100-79784 Lot3 used at a dilution of 1:5,000 (0.2ug/ml) (right). Detection: DAB



Immunoprecipitation: BANP Antibody [NB100-79784] - Detection of human SMAR1 by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-SMAR1 antibody NB100-79784 (lot NB100-79784-3) used for IP at 6 ug per reaction. SMAR1 was also immunoprecipitated by a previous lot of this antibody (lot NB100-79784-2). For blotting immunoprecipitated SMAR1, NB100-79784 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-79784

NBL1-07914	BANP Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-79784

Earn gift cards/discounts by submitting a publication using this product:

www.novusbio.com/publications

