

# Product Datasheet

## ILK Antibody NB100-809

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-809](http://www.novusbio.com/NB100-809)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-809](http://www.novusbio.com/reviews/destination/NB100-809)



**NB100-809**

## ILK Antibody

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA

Product Description	
Description	Novus Biologicals Goat ILK Antibody (NB100-809) is a polyclonal antibody validated for use in IHC, WB and ELISA. Anti-ILK Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Goat
Gene ID	3611
Gene Symbol	ILK
Species	Human, Mouse, Rat, Porcine
Specificity/Sensitivity	All 3 variants (represented by NP_004508.1; NP_001014794.1 and NP_001014795.1) encode the same protein.
Immunogen	Peptide with sequence C-DMIVPILEKMQDK corresponding to C-Terminus according to NP_004508.1, NP_001014794.1, NP_001014795.1.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Peptide ELISA
Recommended Dilutions	Western Blot 0.01-0.03ug/ml, Immunohistochemistry 2-4ug/ml, Immunohistochemistry-Paraffin 2-4ug/ml, Peptide ELISA Detection limit 1:64000
Application Notes	Western blot: Approx. 37kDa band observed in Human, Mouse, Rat and Pig Liver lysate (calculated MW of 34.7kDa according to NP_000036.2). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. IHC: In paraffin embedded Human Kidney shows preferential staining of PCT.

## Images


Western Blot: ILK Antibody [NB100-809] - (0.01ug/ml) staining of Human Liver lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.

250kDa  
150kDa  
100kDa  
75kDa  
50kDa  
37kDa  
25kDa  
20kDa  
15kDa

Immunohistochemistry-Paraffin: ILK Antibody [NB100-809] - (4ug/ml) staining of paraffin embedded Human Kidney. Microwaved antigen retrieval with Tris/EDTA buffer pH9, HRP-staining.



Western Blot: ILK Antibody [NB100-809] - (0.03ug/ml) staining of Pig (A), Mouse (B) and Rat (C) Liver lysates (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.



## Publications

Hannigan GE, Leung-Hagesteijn C, Fitz-Gibbon L et al. Regulation of cell adhesion and anchorage-dependent growth by a new beta 1-integrin-linked protein kinase. Nature 1996-01-04 [PMID: 8538749]



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-809**

---

NBL1-11974	ILK Overexpression Lysate
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-809](http://www.novusbio.com/reviews/submit/NB100-809)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)



