

Product Datasheet

eIF2A [p Ser51] Antibody NB100-81896

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-81896

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-81896



NB100-81896

eIF2A [p Ser51] Antibody


Product Information	
Unit Size	0.05 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS, 50% glycerol, 0.5% BSA
Target Molecular Weight	38 kDa

Product Description	
Description	Novus Biologicals Rabbit eIF2A [p Ser51] Antibody (NB100-81896) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-eIF2A Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	83939
Gene Symbol	EIF2A
Species	Human, Mouse, Rat, Porcine, Canine, Fish
Specificity/Sensitivity	Detects endogenous levels of eIF2A only when phosphorylated at serine 51.
Immunogen	The antiserum was produced against synthesized phosphopeptide derived from human eIF2A around the phosphorylation site of serine 51 (E-L-SP-R-R).


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:2000, Immunohistochemistry 1:100-1:300, Immunocytochemistry/ Immunofluorescence 1:50-200, Immunohistochemistry-Paraffin 1:50-1:100

Images


Western Blot: eIF2A [p Ser51] Antibody [NB100-81896] - Analysis of extracts from Jurkat cells untreated or treated with H₂O₂ using eIF2 alpha (Phospho-Ser51) antibody.




Immunohistochemistry-Paraffin: eIF2A [p Ser51] Antibody [NB100-81896] - Paraffin-embedded human lung carcinoma tissue.



Western Blot: eIF2A [p Ser51] Antibody [NB100-81896] - Western blot analysis of lysates from 1) KB, 2) MCF-7 cells, (Green) primary antibody was diluted at 1: 1000,



Immunohistochemistry-Paraffin: eIF2A [p Ser51] Antibody [NB100-81896] - Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue, using NB100-81896(left) or the same antibody preincubated with blocking peptide (right).



Immunocytochemistry/Immunofluorescence: eIF2A [p Ser51] Antibody [NB100-81896] - Staining of methanol-fixed HeLa cells using NB100-81896



Publications

Wang G, Xu J, Zhao J et al. Arf1-mediated lipid metabolism sustains cancer cells and its ablation induces anti-tumor immune responses in mice Nat Commun 2020-01-10 [PMID: 31924786] (IF/IHC, ICC/IF, Mouse)

Liu YQ, Jia MQ, Xie ZH et al. Arrestins contribute to amyloid beta-induced cell death via modulation of autophagy and the alpha7nACh receptor in SH-SY5Y cells. Sci Rep. 2017-06-13 [PMID: 28611418] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NB100-81896

NBL1-10174	eIF2A Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-81896

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



