

Product Datasheet

RelA/NFkB p65 [p Ser276] Antibody NB100-82086

Unit Size: 0.05 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 6

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-82086

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-82086



NB100-82086

RelA/NFκB p65 [p Ser276] Antibody

Product Information	
Unit Size	0.05 ml
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS, 50% glycerol, 0.5% BSA
Target Molecular Weight	65 kDa
Product Description	
Description	Novus Biologicals Rabbit RelA/NFκB p65 [p Ser276] Antibody (NB100-82086) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-RelA/NFκB p65 Antibody: Cited in 6 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5970
Gene Symbol	RELA
Species	Human, Mouse, Rat
Specificity/Sensitivity	Detects endogenous levels of RelA/NFκB p65 only when phosphorylated at serine 276.
Immunogen	Peptide sequence around phosphorylation site of serine 276(R-P-S(p)-D-R) derived from Human RelA/NFκB p65
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:2000, Immunohistochemistry 1:100-1:300, Immunocytochemistry/ Immunofluorescence 1:50-200, Immunoprecipitation 2-5 ug/mg lysate, Immunohistochemistry-Paraffin 1:50-1:100
Application Notes	RelA/NFκB p65 [p Ser276] antibody validated for WB from a verified customer review.




Images


Immunohistochemistry-Paraffin: RelA/NFkB p65 [p Ser276] Antibody [NB100-82086] - Immunohistochemical analysis of paraffin- embedded human lung carcinoma tissue using NB100-82086



Immunocytochemistry/Immunofluorescence: RelA/NFkB p65 [p Ser276] Antibody [NB100-82086] - Staining of methanol-fixed Hela cells.




Immunohistochemistry-Paraffin: RelA/NFkB p65 [p Ser276] Antibody [NB100-82086] - Human breast carcinoma tissue.




Western Blot: RelA/NFkB p65 [p Ser276] Antibody [NB100-82086] - Western Blot analysis of various cells using NB100-82086 diluted at 1:1000




Immunohistochemistry-Paraffin: RelA/NFkB p65 [p Ser276] Antibody [NB100-82086] - Analysis of paraffin- embedded human breast carcinoma tissue using NB100-82086 incubated with blocking peptide (right)



Immunocytochemistry/Immunofluorescence: RelA/NFkB p65 [p Ser276] Antibody [NB100-82086] - Immunofluorescence staining of methanol-fixed MEF cells untreated or treated with TNF using NB100-82086



Western Blot: RelA/NFkB p65 [p Ser276] Antibody [NB100-82086] - Lane 1: HEK293. Lane 2: Lentivirus mediated NFkB overexpression. We image submitted by a verified customer review.



Publications

Masayuki Harada, Kanae Su-Harada, Takeshi Kimura, Koh Ono, Noboru Ashida Sustained activation of NF- κ B through constitutively active IKK β leads to senescence bypass in murine dermal fibroblasts *Cell Cycle* 2024-03-10 [PMID: 38461418]

Daehong Kim, Giljun Park, Jani Huuhtanen, Bishwa Ghimire, Hanna Rajala, Richard Moriggl, Wing C. Chan, Matti Kankainen, Mikko Myllymäki, Satu Mustjoki STAT3 activation in large granular lymphocyte leukemia is associated with cytokine signaling and DNA hypermethylation *Leukemia* 2021-06-01 [PMID: 34075200]

Gao AY, Diaz Espinosa AM, Gianl F et al. Pim-1 Kinase is a Positive Feedback Regulator of the Senescent Lung Fibroblast Inflammatory Secretome *American journal of physiology. Lung cellular and molecular physiology* 2022-10-12 [PMID: 36223640] (WB, Human)

Kaya S Protective effects of hesperetin on lipopolysaccharide-induced acute lung injury in a rat model *Turk Gogus Kalp Dama* 2020-04-20 [PMID: 30807747] (IF/IHC, IHC-P, Rat)

Zhu Shijun, Moore Terry W, Morii Nao et al. Synthetic curcumin analog UBS109 inhibits the growth of head and neck squamous cell carcinoma xenografts. *Curr Cancer Drug Targets*. 2014-05-01 [PMID: 24628271]

Zhu S, Moore TW, Lin X et al. Synthetic curcumin analog EF31 inhibits the growth of head and neck squamous cell carcinoma xenografts. *Integr Biol (Camb)*. 2012-06-01 [PMID: 22532032] (IHC-P, WB, Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-82086

NBL1-15276	RelA/NFkB p65 Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-82086

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

