

Product Datasheet

HERV Antibody - BSA Free NB100-93579

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 **Publications: 4**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-93579

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-93579



NB100-93579


HERV Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS
Target Molecular Weight	59.9 kDa
Product Description	
Description	Novus Biologicals Rabbit HERV Antibody - BSA Free (NB100-93579) is a polyclonal antibody validated for use in IHC and WB. Anti-HERV Antibody: Cited in 4 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	30816
Gene Symbol	ERVW-1
Species	Human
Specificity/Sensitivity	Human Endogenous Retroviruses (HERVs), a class of multiple-sclerosis-associated retroviruses
Immunogen	Synthetic peptide corresponding to residues 42-50 [QRPGNIDAPC] of the HERV protein.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunohistochemistry-Frozen, In vitro assay, Block/Neutralize
Recommended Dilutions	Western Blot 1:250-1:1000, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, In vitro assay, Block/Neutralize
Application Notes	This antibody is useful for Western Blot. Use in Immunohistochemistry-Paraffin reported in scientific literature (PMID 24635849). Use in In vitro, Immunohistochemistry-Frozen and blocking/neutralizing reported in scientific literature (PMID 26154387) The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.



Images

Western Blot: HERV Antibody [NB100-93579] - WB detection of HERV in human placenta lysate with NB100-93579 antibody at a dilution of 1:250



Publications

Pal D, Patel M, Boulet F et al. H4K16ac activates the transcription of transposable elements and contributes to their cis-regulatory function *Nature Structural & Molecular Biology* 2023-07-01 [PMID: 37308596] (Immunohistochemistry-Paraffin, Human)

Colson A, Depoix CL, Lambert I et al. Specific HIF-2? (Hypoxia-Inducible Factor-2) Inhibitor PT2385 Mitigates Placental Dysfunction In Vitro and in a Rat Model of Preeclampsia (RUPP) Hypertension (Dallas, Tex. : 1979) 2023-05-01 [PMID: 36876500] (WB, Human)

Frese S, Ruebner M, Suhr F et al. Long-Term Endurance Exercise in Humans Stimulates Cell Fusion of Myoblasts along with Fusogenic Endogenous Retroviral Genes In Vivo *PLoS ONE*. 2015-07-09 [PMID: 26154387] (B/N, IHC-Fr, In vitro, Human)

Details:

HERV / Syncytin-1 antibody was used at 1:200 dilution for IHC staining of human muscle cryosections and normal human placental control tissues (third trimester) were employed as control for the staining evaluation (Fig 5). Moreover, the antibody was used for in-vitro blocking/neutralization experiments involving cultured primary human myoblasts (0.5 - 1.0 ug/ml concentration; data in Fig 7, S1).

Buslei R, Strissel PL, Henke C et al. Activation and regulation of endogenous retroviral genes in the human pituitary gland and related endocrine tumors. *Neuropathol. Appl. Neurobiol.* 2014-03-17 [PMID: 24635849] (IHC-P, Human)

Details:

IHC: Normal pituitary adenohypophysis and pituitary adenomas, Figs 1,2; WB: primary cultures of pituitary adenomas, Fig 5A.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-93579

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-93579

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

