

# Product Datasheet

## ARAP1 Antibody (ARAP1-2) [DyLight 488] NB110-68801G

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB110-68801G](http://www.novusbio.com/NB110-68801G)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB110-68801G](http://www.novusbio.com/reviews/destination/NB110-68801G)



**NB110-68801G**

ARAP1 Antibody (ARAP1-2) [DyLight 488]

| Product Information |   |
|---------------------|---|
| Unit Size           | 0.1 ml  |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.  |
| Clonality           | Monoclonal  |
| Clone               | ARAP1-2   |
| Preservative        | 0.05% Sodium Azide  |
| Isotype             | IgG1  |
| Conjugate           | DyLight 488   |
| Purity              | Protein A purified  |
| Buffer              | 50mM Sodium Borate  |

| Product Description     |   |
|-------------------------|---|
| Host                    | Mouse   |
| Gene ID                 | 116985  |
| Gene Symbol             | ARAP1   |
| Species                 | Human   |
| Specificity/Sensitivity | The ARAP1-2 reacts with C-terminal part of ARAP1, a 160 kDa adaptor protein.      |
| Immunogen               | His6-tagged protein corresponding to amino acids 1190-1450 of human ARAP1         |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB110-68801G**

---

|                    |  |
|--------------------|--|
| NBP1-97005G        | Mouse IgG1 Isotype Control (MG1) [DyLight 488] |
| NB110-68801UV      | ARAP1 Antibody (ARAP1-2) [DyLight 350]         |
| H00116985-P01-10ug | Recombinant Human ARAP1 GST (N-Term) Protein   |
| 236-EG-200         | EGF [Unconjugated]                             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB110-68801G](http://www.novusbio.com/reviews/submit/NB110-68801G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

