

Product Datasheet

TfR (Transferrin R) Antibody (23D10) - BSA Free NB110-8468

Unit Size: 0.2 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-8468

Updated 1/22/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-8468



NB110-8468


TfR (Transferrin R) Antibody (23D10) - BSA Free

Product Information	
Unit Size	0.2 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	23D10
Preservative	0.09% Sodium Azide
Isotype	IgG2b
Purity	Protein G purified
Buffer	PBS (pH 7.4)
Product Description	
Description	Novus Biologicals Mouse TfR (Transferrin R) Antibody (23D10) - BSA Free (NB110-8468) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	7037
Gene Symbol	TFRC
Species	Human
Marker	Recycling Endosome Marker
Specificity/Sensitivity	Human soluble transferrin receptor and complex between soluble transferrin receptor and transferrin.
Immunogen	Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with purified human soluble transferrin receptor.
Product Application Details	
Applications	Western Blot, ELISA, Sandwich ELISA
Recommended Dilutions	Western Blot 1:100-1:2000, ELISA 1:100-1:2000, Sandwich ELISA 1:100-1:2000
Application Notes	Recommended pairs for ELISA (capture - detection): NB110-8468 - NB100-73092




Images

ELISA: Transferrin Receptor Antibody (23D10) [NB110-8468] - Titration curve of purified soluble transferrin receptor (sTfR) using pair 23D10-13E4. Capture: MAb 23D10; 200 ng/well Detection: MAb 13E4 conjugated with HRP Room temperature



ELISA: Transferrin Receptor Antibody (23D10) [NB110-8468] - Titration curve of sTfR-Tf complex contained in normal human serum using sTfR-Tf assay. Capture: MAb 23D10; 200 ng/well Detection: MAb 8B9 or 11D3 conjugated with HRP Room temperature





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-8468

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP2-27231	Mouse IgG2b Isotype Control (MPC-11)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-8468

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

