

Product Datasheet

SIAH1/2 Antibody (8G7H12) - BSA Free NB110-93596

Unit Size: 0.5 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-93596

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-93596



NB110-93596

SIAH1/2 Antibody (8G7H12) - BSA Free


Product Information	
Unit Size	0.5 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	8G7H12
Preservative	0.1% Sodium Azide
Isotype	IgG1
Purity	Tissue culture supernatant
Buffer	Tissue culture supernatant
Target Molecular Weight	37 kDa

Product Description	
Description	Novus Biologicals Mouse SIAH1/2 Antibody (8G7H12) - BSA Free (NB110-93596) is a monoclonal antibody validated for use in WB and IP. Anti-SIAH1/2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	6477
Gene Symbol	SIAH1
Species	Human, Rat, Porcine, Drosophila, Zebrafish, Mouse (Negative)
Reactivity Notes	Does not react with mouse.
Specificity/Sensitivity	This is specific for human SIAH 1 and 2 and Drosophila SINA in proliferating cells and stem cells. This does not detect Drosophila SIAH or normal, nondividing cells.
Immunogen	Synthetic peptide made to a C-terminal portion of the drosophila SINA protein (within residues 280-331). [Swiss-Prot# P21461]

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500, Immunoprecipitation 1:10-1:500
Application Notes	<p>This SIAH1/2 antibody is useful for Immunoprecipitation and Western blot, where a band is seen at ~37 kDa.</p> <p>The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.</p>

Images

Western Blot: SIAH1/2 Antibody (8G7H12) [NB110-93596] - Anti-SINA/SIAH monoclonal antibodies recognize both *Drosophila* SINA and human SIAH. There are two SINA-like E3 ligases, SINA and SINAH, in *Drosophila*. SINA monoclonal antibody 8G7H12 is highly speci



Publications

Schmidt, RL et al. Inhibition of RAS-mediated transformation and tumorigenesis by targeting the downstream E3 ubiquitin ligase seven in absentia homologue. *Cancer Res*;67(24):11798-810. 2007-12-15 [PMID: 18089810] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-93596

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NB100-105	HIF-1 alpha Antibody (H1alpha67)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-93596

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

