

Product Datasheet

CaMKII alpha/beta [p Thr286, p Thr287] Antibody (22B1) - BSA Free NB110-96869

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 8

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-96869

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-96869



NB110-96869


CaMKII alpha/beta [p Thr286, p Thr287] Antibody (22B1) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	22B1
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS (pH 7.4), 50% Glycerol
Product Description	
Description	Novus Biologicals Mouse CaMKII alpha/beta [p Thr286, p Thr287] Antibody (22B1) - BSA Free (NB110-96869) is a monoclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. Anti-CaMKII alpha/beta Antibody: Cited in 7 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	815
Gene Symbol	CAMK2A
Species	Human, Mouse, Rat
Specificity/Sensitivity	This antibody is specific for a and B subunits of CaMKII only when they are phosphorylated at Thr-286/287 (in B).
Immunogen	Synthetic peptide
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000, ELISA, Immunohistochemistry 1:100, Immunocytochemistry/ Immunofluorescence 1:1000, Immunoprecipitation, Immunohistochemistry-Paraffin 1:10-1:500
Application Notes	1 ug/ml was sufficient for detection of 0.2 ug CamKII by colorimetric immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary. Use in immunoblotting reported in scientific literature (PMID: 27599291).




Images


Western Blot: CaMKII alpha/beta [p Thr286, p Thr287] Antibody (22B1) [NB110-96869] - Western Blot analysis of Mouse Ventricle lysates showing detection of CaMKII alpha/beta protein using Mouse Anti-CaMKII alpha/beta Monoclonal Antibody, Clone 22B1 (NB110-96869). Primary Antibody: Mouse Anti-CaMKII alpha/beta Monoclonal Antibody (NB110-96869) at 1:1000. Analysis of CaMKII alpha/beta and NFAT phosphorylation in ventricles of 14 day old mice over-expressing CaMK.




Immunocytochemistry/Immunofluorescence: CaMKII alpha/beta [p Thr286, p Thr287] Antibody (22B1) [NB110-96869] - Immunocytochemistry/Immunofluorescence analysis using Mouse Anti-CaMKII alpha/beta Monoclonal Antibody, Clone 22B1 (NB110-96869). Tissue: dissociated hippocampal neurons. Species: Rat. Fixation: Cold 4% paraformaldehyde/0.2% glutaraldehyde in 0.1M sodium phosphate buffer. Primary Antibody: Mouse Anti-CaMKII alpha/beta Monoclonal Antibody (NB110-96869) at 1:1000 for 12 hours at 4C. Secondary Antibody: FITC Goat Anti-Mouse IgG (green) at 1:50 for 30 minutes at RT. Magnification: 10X. Courtesy of: Mary Kennedy, Caltech.



Immunohistochemistry: CaMKII alpha/beta [p Thr286, p Thr287] Antibody (22B1) [NB110-96869] - Immunohistochemistry analysis using Mouse Anti-CaMKII alpha/beta Monoclonal Antibody, Clone 22B1 (NB110-96869). Tissue: colon carcinoma. Species: Human. Fixation: Formalin. Primary Antibody: Mouse Anti-CaMKII alpha/beta Monoclonal Antibody (NB110-96869) at 1:5000 for 12 hours at 4C. Secondary Antibody: Biotin Goat Anti-Mouse at 1:2000 for 1 hour at RT. Counterstain: Mayer Hematoxylin (purple/blue) nuclear stain at 200 I for 2 minutes at RT. Magnification: 40x.



Immunohistochemistry: CaMKII alpha/beta [p Thr286, p Thr287] Antibody (22B1) [NB110-96869] - Immunohistochemistry analysis using Mouse Anti-CaMKII alpha/beta Monoclonal Antibody, Clone 22B1 (NB110-96869). Tissue: backskin. Species: Mouse. Fixation: Bouin's Fixative and paraffin-embedded. Primary Antibody: Mouse Anti-CaMKII alpha/beta Monoclonal Antibody (NB110-96869) at 1:100 for 1 hour at RT. Secondary Antibody: FITC Goat Anti-Mouse (green) at 1:50 for 1 hour at RT. Localization: Muscle, hair follicle, epidermis. Backskin obtained from transgenic mice.



Publications

Lee Y, Kang JS, Ham OJ et al. Gut metabolite trimethylamine N-oxide induces aging-associated phenotype of midbrain organoids for the induced pluripotent stem cell-based modeling of late-onset disease *Frontiers in Aging Neuroscience* 2022-08-16 [PMID: 36051303] (Western Blot, Mouse)

LaMoia TE, Hubbard BT, Guerra MT, Nasiri A et Al. Cytosolic calcium regulates hepatic mitochondrial oxidation, intrahepatic lipolysis, and gluconeogenesis via CAMKII activation *Cell Metab* 2024-08-17 [PMID: 39153480]

Ghane MA, Wei W, Yakout DW et al. Arc ubiquitination regulates endoplasmic reticulum-mediated Ca²⁺ release and CaMKII signaling *Frontiers in cellular neuroscience* 2023-03-14 [PMID: 36998269] (Western Blot, Mouse)

Kim Y, Jang E, Lee Y et al. Miniature Transparent Dopamine Sensor Based on Nanosphere Lithography *Advanced Materials Technologies* 2023-04-02 (WB)

Cai B, Kasikara C, Doran AC et al. MerTK signaling in macrophages promotes the synthesis of inflammation resolution mediators by suppressing CaMKII activity. *Sci Signal* 2018-09-25 [PMID: 30254055] (WB, Mouse)

Ersoy BA, Maner-Smith KM, Li Y et al. Thioesterase-mediated control of cellular calcium homeostasis enables hepatic ER stress. *J. Clin. Invest.* 2017-11-20 [PMID: 29202465] (WB, Mouse)

Conlon DM, Thomas T, Fedotova T et al. Inhibition of apolipoprotein B synthesis stimulates endoplasmic reticulum autophagy that prevents steatosis. *J Clin Invest* 2016-09-06 [PMID: 27599291] (IB)

Jiang L, Xu L, Song Y et al. CaMK II/CREB/Wnt/beta-catenin signaling cascade regulates Angiotensin II-induced podocyte injury and albuminuria. *J Biol Chem* 2013-06-26 [PMID: 23803607] (WB, Mouse)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-96869

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-96869

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

