

Product Datasheet

p38 alpha Antibody (9F12) - BSA Free NB110-96875

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-96875

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-96875



NB110-96875


p38 alpha Antibody (9F12) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	9F12
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS, 50% Glycerol
Product Description	
Description	Novus Biologicals Mouse p38 alpha Antibody (9F12) - BSA Free (NB110-96875) is a monoclonal antibody validated for use in IHC, WB, ELISA and IP. Anti-p38 alpha Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	1432
Gene Symbol	MAPK14
Species	Human, Mouse, Rat
Specificity/Sensitivity	Detects approx 38kDa.
Immunogen	Full length recombinant protein expressed in E.coli cells
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry, Immunoprecipitation, Microarray, Proximity Ligation Assay
Recommended Dilutions	Western Blot 1:1000, ELISA, Immunohistochemistry 1:1000, Immunoprecipitation, Immunohistochemistry-Paraffin 1:1000, Proximity Ligation Assay, Microarray
Application Notes	Detects approx 38kDa protein corresponding to p38a MAPK when loaded with 6 ng of purified p38a by chemiluminescent immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary Antibody. Use in Proximity ligation assay reported in scientific literature (PMID: 30987646).



Images

Western Blot: p38 alpha Antibody (9F12) [NB110-96875] - analysis of Human Cell lysates showing detection of p38 MAPK protein using Mouse Anti-p38 MAPK Monoclonal Antibody, Clone 9F12 . Load: 15 ug protein. Block: 1.5% BSA for 30 minutes at RT. Primary Antibody: Mouse Anti-p38 MAPK Monoclonal Antibody at 1:1000 for 2 hours at RT. Secondary Antibody: Sheep Anti-Mouse IgG: HRP for 1 hour at RT.



Immunohistochemistry: p38 alpha Antibody (9F12) [NB110-96875] - Immunohistochemistry analysis using Mouse Anti-p38 alpha Monoclonal Antibody, Clone 9F12 (NB110-96875). Tissue: colon carcinoma. Species: Human. Fixation: Formalin. Primary Antibody: Mouse Anti-p38 alpha Monoclonal Antibody (NB110-96875) at 1:10000 for 12 hours at 4C. Secondary Antibody: Biotin Goat Anti-Mouse at 1:2000 for 1 hour at RT. Counterstain: Mayer Hematoxylin (purple/blue) nuclear stain at 200 l for 2 minutes at RT. Magnification: 40x.



Publications

SchnOder L, Gasparoni G, NordstrOm K et al. Neuronal deficiency of p38 alpha-MAPK ameliorates symptoms and pathology of APP or Tau-transgenic Alzheimer's mouse models *FASEB J.* 2020-05-31 [PMID: 32475008] (Mouse)

Koch, SR;Choi, H;Mace, EH;Stark, RJ; Toll-like receptor 3-mediated inflammation by p38 is enhanced by endothelial nitric oxide synthase knockdown *Cell Commun. Signal* 2019-04-15 [PMID: 30987646] (PLA, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-96875

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-96875

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

