

Product Datasheet

DNMT3A Antibody (64B1446) - BSA Free NB120-13888

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 3 Publications: 126

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB120-13888

Updated 11/11/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB120-13888




NB120-13888

DNMT3A Antibody (64B1446) - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	64B1446
Preservative	0.05% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein G purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Knockout (KO) Validated Mouse DNMT3A Antibody (64B1446) - BSA Free (NB120-13888) is a monoclonal antibody validated for use in IHC, WB, Flow, ICC/IF and ChIP. Anti-DNMT3A Antibody: Cited in 120 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	1788
Gene Symbol	DNMT3A
Species	Human, Mouse, Rat
Reactivity Notes	Use in Rat reported in scientific literature (PMID:34874218) .
Immunogen	This DNMT3A Antibody (64B1446) was raised against bacteria expressed recombinant mouse Dnmt3a. The epitope was found to lie near the C-terminus (a.a. 705-908), see Chen et (2002) for details.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Chromatin Immunoprecipitation, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Chromatin Immunoprecipitation (ChIP), CyTOF-ready, Immunohistochemistry Whole-Mount, Knockdown Validated, Knockout Validated
Recommended Dilutions	Western Blot, Chromatin Immunoprecipitation reported in scientific literature (PMID 24623306), Flow Cytometry 1 ug per million cells, Immunohistochemistry reported in multiple pieces of scientific literature, Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunohistochemistry-Paraffin 5 ug/ml. Use reported in scientific literature (PMID 22134929), Immunohistochemistry-Frozen reported in scientific literature (PMID 22134929), Immunohistochemistry Whole-Mount reported in scientific literature (PMID 26507142), Chromatin Immunoprecipitation (ChIP), CyTOF-ready, Knockout Validated, Knockdown Validated Validated for Knockdown from CiteAb
Application Notes	Western Blot: Detects a band of approximately 120 kDa (predicted molecular weight: 102 kDa). Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM sodium citrate buffer, pH 6.0 for 10-20 min followed by cooling at RT for 20 min.

Images


Western Blot: DNMT3A Antibody (64B1446) [NB120-13888] - Analysis of (A) Dnmt3a transfected 293 cell lysate and (B) untransfected 293 cell lysate using Dnmt3a antibody at 1 ug/mL.




Immunocytochemistry/Immunofluorescence: DNMT3A Antibody (64B1446) [NB120-13888] - Ntera2 cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X TBS + 0.5% Triton-X100. The cells were incubated with anti-DNMT3A (68B1446) [NB120-13888] at a 1:200 dilution overnight at 4C and detected with an anti-mouse Dylight 488 (Green) at a 1:500 dilution. Actin was detected with Phalloidin 568 (Red) at a 1:200 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective.




Immunohistochemistry-Paraffin: DNMT3A Antibody (64B1446) [NB120-13888] - Detection of DNMT3A on human hepatocellular carcinoma tissue section (NBP2-30221) using 1:100 dilution of DNMT3A antibody (clone 64B1446). The antibody generated a specific cytoplasmic staining in all the cancer cells while some of the cells depicted nuclear staining also.




Immunocytochemistry/Immunofluorescence: DNMT3A Antibody (64B1446) [NB120-13888] - Expression of Dnmt3a in the nuclei of post-mitotic neurons in the olfactory epithelium (OE). DNMT3a (red) is co-expressed with neuron-specific tubulin (green) throughout the development of the olfactory epithelium. Data courtesy of A. Jane Roskams, University of British Columbia.




Flow Cytometry: DNMT3A Antibody (64B1446) [NB120-13888] - Intracellular flow cytometric staining of 1×10^6 CHO (A) and HEK-293 (B) cells using Dnmt3a antibody (dark blue). Isotype control shown in orange. An antibody concentration of $1 \mu\text{g}/1 \times 10^6$ cells was used.




Flow Cytometry: DNMT3A Antibody (64B1446) [NB120-13888] - Using the PE direct conjugate, an intracellular stain was performed on NTERA-2 cells with DNMT3A (64B1446) NB120-13888PE (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of $1 \mu\text{g}/\text{mL}$ for 30 minutes at room temperature. Both antibodies were conjugated to phycoerythrin.




Immunocytochemistry/Immunofluorescence: DNMT3A Antibody (64B1446) [NB120-13888] - DNMT1 reveals a SCMC-like cortical localization with aberrant localization in *Nlrp2tm1a/tm1a* derived oocytes. Whole mount immunofluorescence for DNMT3A in unfertilized control oocytes reveals a characteristic metaphase associated localization and no difference is noted in oocytes derived from *Nlrp2tm1a/tm1a* adams. Image collected and cropped by CiteAb from the following publication (<https://www.nature.com/articles/srep44667>), licensed under a CC-BY license.




Flow Cytometry: DNMT3A Antibody (64B1446) [NB120-13888] - An intracellular stain was performed on Jurkat cells with DNMT3A (64B1446) NB120-13888AF647 (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of $2.5 \mu\text{g}/\text{mL}$ for 30 minutes at room temperature. Both antibodies were conjugated to Alexa Fluor 647.




Western Blot: DNMT3A Antibody (64B1446) [NB120-13888] - NTERA-2 cells were treated with Zebularine as indicated for 24 hours. Cell lysates were prepared and separated on a 7.5% gel by SDS-PAGE. Protein was transferred to PVDF membrane and blocked in 5% non-fat milk. The membrane was probed with 2 ug/mL anti-Dnmt3a in 1% milk, and detected with an anti-mouse HRP secondary antibody using chemiluminescence. Note the decrease in Dnmt3a expression upon treatment with 500 uM Zebularine (arrowhead). Additional bands at 90 and 75 kDa can also be detected with this antibody and may represent alternative splice variants. Tubulin is shown as a loading control.




Flow Cytometry: DNMT3A Antibody (64B1446) [NB120-13888] - An intracellular stain was performed on HepG2 cells with DNMT3A (64B1446) NB120-13888APC (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 1 ug/mL for 30 minutes at room temperature. Both antibodies were conjugated to allophycocyanin (APC).




Flow Cytometry: DNMT3A Antibody (64B1446) [NB120-13888] - An intracellular stain was performed on HepG2 cells with DNMT3A (64B1446) NB120-13888AF647 (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 2.5 ug/mL for 30 minutes at room temperature. Both antibodies were conjugated to Alexa Fluor 647.




Flow Cytometry: DNMT3A Antibody (64B1446) [NB120-13888] - An intracellular stain was performed on HeLa cells with DNMT3A (64B1446) NB120-13888APC (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 1 ug/mL for 30 minutes at room temperature. Both antibodies were conjugated to allophycocyanin (APC).




Western Blot: DNMT3A Antibody (64B1446) - BSA Free [NB120-13888] - DNMT3A knockdown increases miR-200b expression via promoter demethylation. Western blot assay of DNMT3A expression in MDA-MB-231 and BT549 cells after transfection with DNMT3A siRNA (si-DNMT3A) or negative control (NC). Image collected and cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/30324719/>) licensed under a CC-BY license.



Flow Cytometry: DNMT3A Antibody (64B1446) - BSA Free [NB120-13888] - Protein levels of Dnmt3a & MeCP2 are reduced in inhibitory neurons of conditional knockout mice. (A) Western blot of half brain hemisphere from *Dnmt3a* or *Mecp2* cKO mice demonstrating loss of each protein. (B) Immunofluorescence (IF) images of WT, *Dnmt3a* cKO or *Mecp2* cKO mice probing for *Dnmt3a* or MeCP2 (red), the Sun1-sfGFP-myc-PA fusion protein that marks the nuclear envelope & is dependent on Cre expression (green), & DAPI to mark genomic DNA (blue). $n = 3$ mice per genotype (western), $n = 3$ mice (IF) representative images for two brain regions shown. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32159514/>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



DNMT1 reveals a SCMC-like cortical localization with aberrant localization in *Nlrp2tm1a/tm1a* derived oocytes. (A) Whole mount immunofluorescence co-staining with anti-DNMT1 and DAPI on the same *Nlrp2M+/Z+* embryos used to stain NLRP2 reveals a SCMC-like cortical localization of DNMT1 with a diffuse nuclear stain in later stage embryos. Scale bar represented is 200 μ m. (B) DNMT1 staining in paraffin embedded oocyte sections reveals a less cortical, more cytoplasmic intense stain with nucleolar focus in oocytes from *Nlrp2tm1a/tm1a* dams compared to oocytes from *Nlrp2+/+* dams. (C) Whole mount immunofluorescence for DNMT3A in unfertilized control oocytes reveals a characteristic metaphase associated localization and no difference is noted in oocytes derived from *Nlrp2tm1a/tm1a* dams. (D) Whole mount immunofluorescence in unfertilized oocytes reveals that as expected, DNMT3B is not yet expressed in unfertilized control oocytes and no expression is noted in oocytes derived from *Nlrp2tm1a/tm1a* dams. Image collected and cropped by CiteAb from the following open publication (<https://pubmed.ncbi.nlm.nih.gov/28317850/>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Chatterjee B, Majumder P, Chen C et al. Hypoxia-induced genome-wide DNA demethylation by DNMT3A and EMT of cancer cells Cellular & Molecular Biology Letters 2025-08-05 [PMID: 40764968]

Kong X, Chen J, Xie W, Brown SM et Al. Defining UHRF1 Domains that Support Maintenance of Human Colon Cancer DNA Methylation and Oncogenic Properties Cancer Cell 2019-04-09 [PMID: 30956060]

Masaki Yagi, Satoshi Kishigami, Akito Tanaka, Katsunori Semi, Eiji Mizutani, Sayaka Wakayama, Teruhiko Wakayama, Takuya Yamamoto, Yasuhiro Yamada Derivation of ground-state female ES cells maintaining gamete-derived DNA methylation. Nature 2017-12-14 [PMID: 28746308]

Nguyet-Minh Hoang, Yunxia Liu, Paul D. Bates, Alexa R. Heaton, Angelica F. Lopez, Peng Liu, Fen Zhu, Ruoyu Chen, Apoorv Kondapelli, Xiyu Zhang, Paul E. Selberg, Vu N. Ngo, Melissa C. Skala, Christian M. Capitini, Lixin Rui Targeting DNMT3A-mediated oxidative phosphorylation to overcome ibrutinib resistance in mantle cell lymphoma Cell Reports Medicine 2024-03-29 [PMID: 38554704]

Lisa-Marie Brenner, Florian Meyer, Haiqian Yang, Anja R Köhler, Pavel Bashtrykov, Ming Guo, Albert Jeltsch, Cristiana Lungu, Monilola A Olayioye Repeat DNA methylation is modulated by adherens junction signaling. Communications biology 2024-03-11 [PMID: 38454140]

Yang C, Deng L, Bao F et al. Sevoflurane with Low Concentration Decrease DNA Methylation on Chronic Obstructive Pulmonary Disease (COPD)-Related Gene Promoter in COPD Rat COPD 2023-12-01 [PMID: 38010369] (WB)

Kubo N, Uehara R, Uemura S et al. Combined and differential roles of ADD domains of DNMT3A and DNMT3L on DNA methylation landscapes in mouse germ cells bioRxiv 2023-09-08 (ICC/IF, Mouse)

Ren G, Li H, Hong D et al. LINC00955 suppresses colorectal cancer growth by acting as a molecular scaffold of TRIM25 and Sp1 to Inhibit DNMT3B-mediated methylation of the PHIP promoter BMC cancer 2023-09-23 [PMID: 37742010] (WB, Human)

Narabayashi H, Koma C, Nakata K et al. Gut microbiota-dependent adaptor molecule recruits DNA methyltransferase to the TLR4 gene in colonic epithelial cells to suppress inflammatory reactions Frontiers in molecular biosciences 2022-10-21 [PMID: 36339704] (ChIP, Mouse)

Yu V, Yong F, Chen K et al. Establishment of beta cell heterogeneity via differential CpG methylation atNnat bioRxiv 2023-02-05 (IHC, Mouse)

da C. Fernandes C, da Silva R, de Almeida G et al. Epigenetic Differences Arise in Endothelial Cells Responding to Cobalt-Chromium Journal of Functional Biomaterials 2023-02-26 [PMID: 36976051] (WB, Human)

Fernandes C, da Silva R, Wood P et al. Titanium-Enriched Medium Promotes Environment-Induced Epigenetic Machinery Changes in Human Endothelial Cells Journal of Functional Biomaterials 2023-02-27 [PMID: 36976055]

More publications at <http://www.novusbio.com/NB120-13888>



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB120-13888

NBL1-09970	DNMT3A Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB120-13888

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

