

Product Datasheet

n-Myc Antibody (NMYC-1) - BSA Free NB200-109

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 6

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB200-109

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB200-109



NB200-109


n-Myc Antibody (NMYC-1) - BSA Free

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | 1.0 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | NMYC-1 |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG2a Kappa |
| Purity | Protein A purified |
| Buffer | PBS |
| Product Description | |
| Description | Novus Biologicals Mouse n-Myc Antibody (NMYC-1) - BSA Free (NB200-109) is a monoclonal antibody validated for use in WB, Flow, ICC/IF, IP and ChIP. Anti-n-Myc Antibody: Cited in 6 publications. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Mouse |
| Gene ID | 4613 |
| Gene Symbol | MYCN |
| Species | Human, Mouse |
| Specificity/Sensitivity | This antibody does not react with c-Myc or L-Myc. |
| Immunogen | A synthetic peptide (SPYVESEDAPPQKC) conjugated to KLH, corresponding to amino acids 327-339 of human n-Myc. [UniProt# P04198] |
| Product Application Details | |
| Applications | Western Blot, Chromatin Immunoprecipitation, Flow Cytometry, Gel Super Shift Assays, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), CyTOF-ready |
| Recommended Dilutions | Western Blot 1:500 - 1:1000, Chromatin Immunoprecipitation 1:10 - 1:500. Use reported in scientific literature (PMID 21068393), Flow Cytometry reported in scientific literature (PMID 27539729), Immunocytochemistry/ Immunofluorescence 1:50 - 1:100, Immunoprecipitation 1:10 - 1:100, Gel Super Shift Assays 1:1 - 1:100. Use reported in scientific literature (PMID 21068393), Chromatin Immunoprecipitation (ChIP) 1:10-1:500, CyTOF-ready |
| Application Notes | In Western blot a band can be seen at ~50-65 kDa. |




Images


Western Blot: n-Myc Antibody (NMYC-1) [NB200-109] - Analysis of n-Myc expression in IMR32 whole cell lysate.



Immunocytochemistry/Immunofluorescence: n-Myc Antibody (NMYC-1) [NB200-109] - The n-Myc (NMYC-1) antibody was tested in NTERA2 cells at a 1:50 dilution against Dylight 488 (Green). Actin and nuclei were counterstained against Phalloidin 568 (Red) and DAPI (Blue), respectively.



Western Blot: n-Myc Antibody (NMYC-1) [NB200-109] - Analysis of n-Myc expression in human heart tissue.



Publications

Ton AT, Foo J, Singh K et al. Development of VPC-70619, a Small-Molecule N-Myc Inhibitor as a Potential Therapy for Neuroendocrine Prostate Cancer International journal of molecular sciences 2022-02-26 [PMID: 35269731] (PLA, Human)

Gupta R, Turati V, Brian D et al. Nov/CCN3 Enhances Cord Blood Engraftment by Rapidly Recruiting Latent Human Stem Cell Activity Cell Stem Cell. [PMID: 32197066] (FLOW, Human)

Yang Y, Ding L, Zhou Q et al. Silencing of AURKA augments the antitumor efficacy of the AURKA inhibitor MLN8237 on neuroblastoma cells Cancer Cell Int. 2020-01-07 [PMID: 31920463] (WB, Mouse)

Murakami-Tonami Y, Ikeda H, Yamagishi R et al. SGO1 is involved in the DNA damage response in MYCN-amplified neuroblastoma cells. Sci Rep. 2016-08-19 [PMID: 27539729] (WB, FLOW, Chemotaxis, Human)

Murakami-Tonami Y, Kishida S, Takeuchi I et al. Inactivation of SMC2 shows a synergistic lethal response in MYCN-amplified neuroblastoma cells. Cell Cycle 2014-01-01 [PMID: 24553121]

Tasseva G, Cole L, Vance JE. N-Myc and SP regulate phosphatidylserine synthase-1 expression in brain and glial cells. J Biol Chem 286(2):1061-73. 2011-01-14 [PMID: 21068393]



Procedures

Western Blot Protocol for NB200-109: n-Myc Antibody (NMYC-1)

n-Myc Antibody (NMYC-1):

Western Blot Protocol

1. Perform SDS-PAGE on samples to be analyzed, loading 40 ug of total protein per lane.
2. Transfer proteins to membrane according to the instructions provided by the manufacturer of the membrane and transfer apparatus.
3. Stain according to standard Ponceau S procedure (or similar product) to assess transfer success, and mark molecular weight standards where appropriate.
4. Rinse the blot.
5. Block the membrane using standard blocking buffer for at least 1 hour.
6. Wash the membrane in wash buffer three times for 10 minutes each.
7. Dilute primary antibody in blocking buffer and incubate 1 hour at room temperature.
8. Wash the membrane in wash buffer three times for 10 minutes each.
9. Apply the diluted HRP conjugated secondary antibody in blocking buffer (as per manufacturers instructions) and incubate 1 hour at room temperature.
10. Wash the blot in wash buffer three times for 10 minutes each (this step can be repeated as required to reduce background).
11. Apply the detection reagent of choice in accordance with the manufacturers instructions.

Note: Tween-20 can be added to the blocking or antibody dilution buffer at a final concentration of 0.05-0.2%.

*The above information is only intended as a guide. The researcher should determine what protocol best meets their needs. Please follow safe laboratory procedures.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB200-109

| | |
|------------------|--|
| NBL1-13417 | n-Myc Overexpression Lysate |
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP] |
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB7539 | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-96981-0.5mg | Mouse IgG2a Kappa Isotype Control (M2AK) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB200-109

Earn gift cards/discounts by submitting a publication using this product:

www.novusbio.com/publications

