

Product Datasheet

DDX3 Antibody - BSA Free NB200-196

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB200-196

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB200-196



NB200-196

DDX3 Antibody - BSA Free


Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)
Target Molecular Weight	73 kDa

Product Description	
Description	Novus Biologicals Rabbit DDX3 Antibody - BSA Free (NB200-196) is a polyclonal antibody validated for use in WB and IP. Anti-DDX3 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1654
Gene Symbol	DDX3X
Species	Human
Reactivity Notes	Orangutan, Chimpanzee (100%).
Immunogen	The immunogen recognized by this antibody maps to a region between residue 600 and the C-terminus (residue 661) of human DEAD (Asp-Glu-Ala-Asp) box polypeptide 3 using the numbering given in entry NP_001347.2 (GeneID 1654).

Product Application Details	
Applications	Western Blot, Immunoprecipitation, Knockdown Validated
Recommended Dilutions	Western Blot 1:5000-1:15000, Immunoprecipitation 2-5 ug/mg lysate, Knockdown Validated

Images

Western Blot: DDX3 Antibody [NB200-196] - Detection of DDX3 by Western Blot and Immunoprecipitation. Samples: A. Whole cell lysate from HeLa (5, 15 and 50 mcg) and Ramos (50 mcg) cells. B. Whole cell lysate (1 mg/IP; 1/4 of IP loaded/lane) from HeLa cells. Antibodies: NB 200-196 used at 0.04 mcg/ml (A) and 0.1 mcg/ml (B) for Western Blot and at 3 mcg/mg lysate for IP (B). DDX3 was also immunoprecipitated using NB 200-194 and NB 200-195 at 3 mcg/mg lysate. Detection: Chemiluminescence with an exposure time of 10 seconds (A and B).



Publications

Saikruang, W;Ang Yan Ping, L;Abe, H;Kasumba, DM;Kato, H;Fujita, T; The RNA helicase DDX3 promotes IFNB transcription via enhancing IRF-3/p300 holocomplex binding to the IFNB promoter Scientific reports [PMID: 35273248] (Mouse)

Nguyen CN, Nguyen KVA, Eom S et al. Regulation of p21 expression for anti-apoptotic activity of DDX3 against sanguinarine-induced cell death on intrinsic pathway Phytomedicine 2019-09-23 [PMID: 31568920] (WB, KD, Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB200-196

NB800-PC1	HeLa Whole Cell Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB200-196

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

