

Product Datasheet

GAP-43 [p Ser41] Antibody - Azide Free NB300-179

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-179

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-179




NB300-179**GAP-43 [p Ser41] Antibody - Azide Free**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	10mM HEPES (pH 7.5), 0.15M NaCl, 0.1 mg/ml BSA and 50% Glycerol
Target Molecular Weight	50 kDa
Product Description	
Description	Novus Biologicals Rabbit GAP-43 [p Ser41] Antibody - Azide Free (NB300-179) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	2596
Gene Symbol	GAP43
Species	Mouse, Rat
Marker	Neuronal Marker
Specificity/Sensitivity	Specific for endogenous levels of the ~50 kDa Gap-43 protein phosphorylated at Ser41. In some tissues the antibody also recognizes a higher molecular weight protein that is also recognized by the pan Gap-43 antibody, which may be a Gap-43 aggregate or oligomer. Immunolabeling is completely eliminated by treatment with lambda-phosphatase.
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser41 conjugated to KLH. Accession # P07936
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:1000
Application Notes	NB300-179 can be used in Western blot where a band is seen at ~ 43 kDa, representing phosphorylated Gap43. In some tissues the antibody also recognizes a lower MW band that is also recognized by the PAN Gap43 antibody and, thus, it is likely to be a breakdown product. Immunohistochemical analysis has not yet been determined.

Images

Western Blot: GAP43 [p Ser41] Antibody [NB300-179] - Western blot of rat cortex lysate showing specific immunolabeling of the ~50k Gap-43 protein phosphorylated at Ser41 (Control). The phosphospecificity of this labeling is shown in the second lane (lambda-phosphatase). The blot is identical to the control except that it was incubated in lambda-phosphatase (1200 units for 30 min) before being exposed to the GAP-43 Ser41 antibody. The immunolabeling of GAP-43 is completely eliminated by treatment with lambda-phosphatase.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-179

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-179

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

