

# Product Datasheet

## Dynamin [p Ser778] Antibody - Azide Free NB300-210

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB300-210](http://www.novusbio.com/NB300-210)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB300-210](http://www.novusbio.com/reviews/destination/NB300-210)




**NB300-210****Dynamin [p Ser778] Antibody - Azide Free**

<b>Product Information</b>	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at -20C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	No Preservative
<b>Isotype</b>	IgG
<b>Purity</b>	Antigen Affinity-purified
<b>Buffer</b>	10mM HEPES (pH 7.5), 0.15M NaCl, 0.1 mg/ml BSA and 50% Glycerol
<b>Target Molecular Weight</b>	95 kDa
<b>Product Description</b>	
<b>Description</b>	Novus Biologicals Sheep Dynamin [p Ser778] Antibody - Azide Free (NB300-210) is a polyclonal antibody validated for use in WB. Anti-Dynamin Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Sheep
<b>Gene ID</b>	1759
<b>Gene Symbol</b>	DNM1
<b>Species</b>	Mouse, Rat
<b>Reactivity Notes</b>	Mouse reactivity reported in scientific literature (PMID: 26658127)
<b>Specificity/Sensitivity</b>	Specific for endogenous levels of the ~95 kDa dynamin protein phosphorylated at Ser778. Labels the purified protein phosphorylated in vitro by cdk5 but not by PKC. Does not cross react with other purified substrates of cdk5 (e.g. amphiphysin and synapsin). Immunolabeling is completely eliminated by treatment with lambda-phosphatase.
<b>Immunogen</b>	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser778 conjugated to KLH. Accession # P21575
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot
<b>Recommended Dilutions</b>	Western Blot 1:1000
<b>Application Notes</b>	This antibody is useful in Western blot where a band is seen at ~ 95kD, representing Dynamin. Immunohistochemistry and Immunofluorescence have not been tested.



## Images

Western Blot: Dynamin [p Ser778] Antibody [NB300-210] - Analysis of rat hippocampal lysate stimulated with forskolin showing specific immunolabeling of the ~95k dynamin phosphorylated at Ser778 (Control). The phosphospecificity of this labeling is shown in the second lane (lambda-phosphatase). The blot is identical to the control except that it was incubated in lambda-Ptase (1200 units for 30 min) before being exposed to the Anti-Ser778 dynamin. The immunolabeling is completely eliminated by treatment.



## Publications

Du JJ, Yan L, Zhang W et al. Clathrin-independent but dynamin-dependent mechanisms mediate Ca<sup>2+</sup>-triggered endocytosis of the glutamate GluK2 receptor upon excitotoxicity *Journal of Integrative Neuroscience* 2020-09-30 [PMID: 33070524]

Patel A, Yamashita N, Ascano M et al. RCAN1 links impaired neurotrophin trafficking to aberrant development of the sympathetic nervous system in Down syndrome. *Nat Commun.* 2015-12-15 [PMID: 26658127] (WB, Mouse)

### Details:

phospho-Dynamin [Ser778] antibody was used for WB analysis of lysates from salivary glands of mouse pups P1 (mouse model of Down syndrome) and NGF stimulated sympathetic neurons of P0.5 Sprague–Dawley rats. The lysates were prepared in Laemmli buffer (62.5 mM Tris-HCl, 2% SDS, 10% Glycerol, 5% 2-mercaptoethanol, 0.01% bromophenol blue) and the antibody was employed at 1:5000 dilution.

Tan, TC et al. Cdk5 is Essential for Synaptic Vesicle Endocytosis. *Nat. Cell Bio.*, 5(8): 701-710. 2003-01-01 [PMID: 12855954]



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB300-210**

---

HAF016	Donkey anti-Sheep IgG Secondary Antibody [HRP]
NL010	Donkey anti-Sheep IgG Secondary Antibody [NL557]
NBP1-97055-10mg	Sheep IgG Isotype Control
NBP2-76558PEP	Dynamin Recombinant Protein Antigen

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB300-210](http://www.novusbio.com/reviews/submit/NB300-210)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

