

# Product Datasheet

## IGF-II R/IGF2R Antibody (2G11) - BSA Free NB300-514

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

**Publications: 14**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB300-514](http://www.novusbio.com/NB300-514)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB300-514](http://www.novusbio.com/reviews/destination/NB300-514)



**NB300-514**


IGF-II R/IGF2R Antibody (2G11) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2G11
Preservative	0.05% Sodium Azide
Isotype	IgG2a Kappa
Purity	Protein A purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Mouse IGF-II R/IGF2R Antibody (2G11) - BSA Free (NB300-514) is a monoclonal antibody validated for use in IHC, WB, ELISA, Flow, ICC/IF and IP. Anti-IGF-II R/IGF2R Antibody: Cited in 13 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3482
Gene Symbol	IGF2R
Species	Human, Mouse, Rat, Bovine, Primate
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 9119492). Does not react with Chinese hamster ovary cell Mannose 6 Phosphate Receptor (PMID: 8163521).
Marker	Late Endosome Marker
Immunogen	Purified bovine Mannose 6 Phosphate Receptor (Cation independent) [UniProt# P08169]
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, CyTOF-ready
Recommended Dilutions	Western Blot 1 ug/ml, Flow Cytometry 1 ug per million cells, ELISA, Immunohistochemistry 1:100, Immunocytochemistry/ Immunofluorescence 1:10-1:250, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:100, CyTOF-ready
Application Notes	By Western Blot, a 300 kDa protein representing Mannose 6 Phosphate Receptor is seen in HeLa cell extract under non-reducing conditions. Immunofluorescence staining of Mannose 6 Phosphate Receptor in HeLa cells with this antibody results in perinuclear staining. This antibody is also useful in Flow Cytometry, Immunoprecipitation, ELISA and Immunohistochemistry-on paraffin embedded sections. This antibody is CyTOF ready.




## Images


**Immunocytochemistry/Immunofluorescence: IGF-II R/IGF2R Antibody (2G11) [NB300-514]** - HeLa cells were fixed in 4% paraformaldehyde for 10 minutes and permeabilized in 0.05% Triton X-100 in PBS for 5 minutes. The cells were incubated with anti-IGF-II R/IGF2R Antibody (2G11) NB300-514 at 1 ug/ml overnight at 4C and detected with an anti-mouse DyLight 488 (Green) at a 1:1000 dilution for 60 minutes. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 100X objective and digitally deconvolved.




**Immunohistochemistry: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) [NB300-514]** - Staining of human testis using NB300-514 antibody at 1:100.




**Flow Cytometry: IGF-II R/IGF2R Antibody (2G11) [NB300-514]** - An intracellular stain was performed on MCF7 cells with IGF-II R/IGF2R [2G11] Antibody NB300-514AF647 (blue) and a matched isotype control (orange). Cells were fixed with 4% PFA and then permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 5 ug/mL for 30 minutes at room temperature. Both antibodies were conjugated to Alexa Fluor 647.




**Immunocytochemistry/Immunofluorescence: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) [NB300-514]** - Mannose 6-Phosphate Receptor staining (green), F-Actin staining with Phalloidin (red) and nuclei with DAPI (blue) is shown. Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with or an antibody recognizing Mannose 6-Phosphate Receptor at a dilution of 1:20 over night at 4C, washed with PBS and incubated with a DyLight-488 conjugated.




Immunocytochemistry/Immunofluorescence: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) [NB300-514] - HMVEC Cells




Immunocytochemistry/Immunofluorescence: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) [NB300-514] - Staining in HeLa Cells. Mannose 6-Phosphate Receptor staining (green), F-Actin staining with Phalloidin (red) and nuclei with DAPI (blue) is shown. Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with an antibody recognizing Mannose 6-Phosphate Receptor at a dilution of 1:20 overnight at 4 C, washed with PBS and incubated with a DyLight-488 conjugated secondary antibody.



Immunocytochemistry/Immunofluorescence: IGF-II R/Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) [NB300-514] - Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) [NB300-514] - Mannose 6 Phosphate Receptor (Cation independent) antibody (2G11) was tested at (1:250) in HeLa cells with DyLight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and DyLight-550 (red).



Flow Cytometry: Mannose 6 Phosphate Receptor (Cation independent) Antibody (2G11) [NB300-514] - Intracellular flow cytometric staining of  $1 \times 10^6$  MCF-7 cells using Mannose 6 Phosphate Receptor antibody (dark blue). Isotype control shown in orange. An antibody concentration of  $1 \mu\text{g}/1 \times 10^6$  cells was used.



## Publications

Bird, LE;Xu, B;Hobbs, AD;Ziegler, AR;Scott, NE;Newton, P;Thomas, DR;Edgington-Mitchell, LE;Newton, HJ; Coxiella burnetii manipulates the lysosomal protease cathepsin B to facilitate intracellular success Nature communications 2025-04-24 [PMID: 40274809]

Wang Q, Tao C, Hannan A et al. Lacrimal gland budding requires PI3K-dependent suppression of EGF signaling Science Advances 2021-07-02 [PMID: 34193412]

Calcagni' A, Staiano L, Zampelli N et al. Loss of the batten disease protein CLN3 leads to mis-trafficking of M6PR and defective autophagic-lysosomal reformation Nature communications 2023-07-03 [PMID: 37400440]

Choi WY, Kim S, Aurass P, Huo W et Al. SdhA blocks disruption of the Legionella-containing vacuole by hijacking the OCRL phosphatase Cell Rep 2021-11-03 [PMID: 34731604]

Aguilar-Gonzalez A, Gonzalez-Correa JE, Barriocanal-Casado E et al. Isogenic GAA-KO Murine Muscle Cell Lines Mimicking Severe Pompe Mutations as Preclinical Models for the Screening of Potential Gene Therapy Strategies International journal of molecular sciences 2022-06-04 [PMID: 35682977] (FLOW, Mouse)

Muzio L, Sirtori R, Gornati D et al. Retromer stabilization results in neuroprotection in a model of Amyotrophic Lateral Sclerosis Nat Commun 2020-07-31 [PMID: 32737286] (IF/IHC, Mouse)

Kakegawa W, Katoh A, Narumi S et al. Optogenetic Control of Synaptic AMPA Receptor Endocytosis Reveals Roles of LTD in Motor Learning Neuron 2018-08-08 [PMID: 30122381] (ICC/IF, Mouse)

Takahashi K, Mashima H, Miura K et al. Disruption of Small GTPase Rab7 Exacerbates the Severity of Acute Pancreatitis in Experimental Mouse Models. Sci Rep. 2017-06-06 [PMID: 28588238] (IF/IHC, WB, Mouse)

Barker LP, George KM, Falkow S, Small PL. Differential trafficking of live and dead Mycobacterium marinum organisms in macrophages. Infect Immun. 1997-04-01 [PMID: 9119492] (ICC/IF, Mouse)

Ganley IG, Carroll K, Bittova L, Pfeffer S. Rab9 GTPase regulates late endosome size and requires effector interaction for its stability. Mol Biol Cell. 2004-12-01 [PMID: 15456905] (ICC/IF, WB, Human)

Arighi CN, Hartnell LM, Aguilar RC et al. Role of the mammalian retromer in sorting of the cation-independent mannose 6-phosphate receptor. J Cell Biol. 2004-04-01 [PMID: 15078903] (ICC/IF, Human)

Valdez AC, Cabaniols JP, Brown MJ, Roche PA. Syntaxin 11 is associated with SNAP-23 on late endosomes and the trans-Golgi network. J Cell Sci. 1999-03-01 [PMID: 10036234] (ICC/IF, Human)

More publications at <http://www.novusbio.com/NB300-514>





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB300-514**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB300-514](http://www.novusbio.com/reviews/submit/NB300-514)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

