

Product Datasheet

TGN38 Antibody (2F7.1) NB300-575

Unit Size: 100 uL

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 **Publications: 18**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-575

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-575




NB300-575

TGN38 Antibody (2F7.1)


Product Information	
Unit Size	100 uL
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2F7.1
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein A or G purified
Buffer	PBS, 1 mg/mL BSA
Product Description	
Description	Novus Biologicals Mouse TGN38 Antibody (2F7.1) (NB300-575) is a monoclonal antibody validated for use in IHC, WB, Flow and ICC/IF. Anti-TGN38 Antibody: Cited in 16 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10618
Gene Symbol	TGOLN2
Species	Human, Mouse, Rat, Chinese Hamster, Monkey, Rabbit
Reactivity Notes	Chinese Hamster reactivity reported in scientific literature (PMID: 11408588). Rabbit reactivity reported in scientific literature (PMID: 10514494). Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information. Use in Human reported in scientific literature (PMID:24576880).
Marker	TGN Marker
Specificity/Sensitivity	Detects Trans Golgi Network (TGN) 38 from human and monkey samples and both recombinant and endogenous rat samples.
Immunogen	Synthetic peptide corresponding to residues L(18) P S A S K P N N T S S E N N P P(34) C of rat TGN38.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500, Flow Cytometry 1:20 - 1:50, Immunohistochemistry 1:1000, Immunocytochemistry/ Immunofluorescence 1:1000, Immunohistochemistry-Paraffin 1:1000
Application Notes	IP usage was reported in scientific literature (PMID: 10514494). WB: Detects an approx. 38 kDa protein representing deglycosylated recombinant rat TGN38. Under normal WB conditions, detects an 85-95 kDa band depending on cell type. IF: Staining of TGN38 in NRK cells results in staining predominantly in the TGN. IF: it has been shown that methanol works best as a fixative.

Images


Immunocytochemistry/Immunofluorescence: TGN38 Antibody (2F7.1) [NB300-575] - HeLa and COS7 cells were rinsed with PBS and fixed with -20C methanol for 5 minutes. Cells were incubated with anti-TGN38 antibody diluted 1:100 in PBS containing 0.2% BSA for 1 hour. The cells were then rinsed in PBS and incubated with goat anti-mouse IgG conjugated to rhodamine for 30 minutes in PBS containing 0.2% BSA. Following washing, cells were imaged using a Zeiss Axioplan microscope. ICC/IF image submitted by a verified customer review.




Immunocytochemistry/Immunofluorescence: TGN38 Antibody (2F7.1) [NB300-575] - Staining of TGN-38 in NRK cells.



Immunocytochemistry/Immunofluorescence: TGN38 Antibody (2F7.1) [NB300-575] - Analysis of Phalloidin (green) and TGN38 (purple) in U2OS cells. Formalin fixed cells were permeabilized with 0.1% Triton X-100 in PBS for 10 minutes at room temperature and blocked with 2% BSA in PBS + 0.1% Triton X-100 for 30 minutes at room temperature. Cells were probed with a TGN38 monoclonal antibody at a dilution of 1:75 for at least 1 hour at room temperature, washed with PBS, and incubated with DyLight 680 goat anti-mouse IgG secondary antibody at a dilution of 1:250 for 30 minutes at room temperature. Actin was stained with DyLight 488 Phalloidin at a dilution of 1:300 (1 unit/mL final concentration) for 30 minutes. Images were taken on a Thermo Scientific ArrayScan VTI at 20X magnification.



Immunocytochemistry/Immunofluorescence: TGN38 Antibody (2F7.1) [NB300-575] - Staining in NS-1 Cells.



Publications

Koepke J, Long W, Barr A et al. Furin may contribute to proglucagon processing and glucagon-like Peptide-1 production in human alpha cells. *Molecular metabolism* 2025-09-27 [PMID: 41022304]

Boyer CK, Bauchle CJ, Zhang J et al. Synchronized proinsulin trafficking reveals delayed Golgi export accompanies β -cell secretory dysfunction in rodent models of hyperglycemia *Scientific reports* 2023-03-30 [PMID: 36997560] (ICC/IF, Rat)

Parchure A, Tian M, Stalder D et al. Liquid-liquid phase separation facilitates the biogenesis of secretory storage granules *The Journal of cell biology* 2022-12-05 [PMID: 36173346] (Immunocytochemistry/ Immunofluorescence, Rat)

Andreska T, Lüningschrör P, Wolf D et al. DRD1 signaling modulates TrkB turnover and BDNF sensitivity in direct pathway striatal medium spiny neurons *Cell reports* 2023-05-29 [PMID: 37252844] (IHC, Rat)

J Tisdale E, R Artalejo C Rab2 stimulates LC3 lipidation on secretory membranes by noncanonical autophagy *Experimental cell research* 2023-05-16 [PMID: 37201743]

Maruo K, Nishiyama M, Honda Y et al. Increased GLUT1 expression and localization to Golgi apparatus of acinar cells in the parotid gland of Goto-Kakizaki diabetic rats *Archives of oral biology* 2022-12-05 [PMID: 36521282] (IHC-P, Rat)

Rohli KE, Boyer CK, Bearrows SC et al. ER Redox Homeostasis Regulates Proinsulin Trafficking and Insulin Granule Formation in the Pancreatic Islet beta-Cell Function (Oxford, England) 2022-09-28 [PMID: 36325514] (IF/IHC, Mouse)

Kamat, V;Radtke, JR;Hu, Q;Wang, W;Sweet, IR;Hampe, CS; Autoantibodies directed against glutamate decarboxylase interfere with glucose-stimulated insulin secretion in dispersed rat islets *International journal of experimental pathology* [PMID: 35246889]

Andreska T Effects of dopamine on BDNF/TrkB mediated signaling and plasticity on cortico-striatal synapses Thesis 2021-01-01 (IF/IHC, Rat)

Li M, Shao F, Qian Q Et al. A putative long noncoding RNA-encoded micropeptide maintains cellular homeostasis in pancreatic beta cells *Molecular Therapy - Nucleic Acids* 2021-12-01 [PMID: 34513312] (ICC/IF, Rat)

Alsalem J A, Patel D et al. Characterization of vitamin D production by human ocular barrier cells. *Invest Ophthalmol Vis Sci* 2014-07-04 [PMID: 24576880] (ICC/IF, Human)

Casas M, FadO R, DomIniguez JL et al. Sensing of nutrients by CPT1C controls SAC1 activity to regulate AMPA receptor trafficking *J. Cell Biol.* 2020-10-05 [PMID: 32931550] (ICC/IF, Mouse)

More publications at <http://www.novusbio.com/NB300-575>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-575

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-575

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

