

Product Datasheet

Leptin/OB Antibody NB300-611

Unit Size: 50 uL

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-611

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-611



NB300-611


Leptin/OB Antibody

Product Information	
Unit Size	50 uL
Concentration	2 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS with 1 mg/ml BSA
Product Description	
Description	Novus Biologicals Rabbit Leptin/OB Antibody (NB300-611) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Leptin/OB Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3952
Gene Symbol	LEP
Species	Human, Mouse, Rat, Porcine, Avian, Sheep
Reactivity Notes	Rat reactivity reported in scientific literature (PMID: 23270329). Human reactivity reported in scientific literature (PMID: 21983012). Porcine reactivity reported in scientific literature (PMID: 18413648). Avian reactivity reported in scientific literature (PMID: 16730725).
Specificity/Sensitivity	Leptin
Immunogen	Synthetic peptide corresponding to residues Q(25) K V Q D D T K T L I K T I V T R I N D(44) of mouse Leptin.
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Frozen, Immunocytochemistry
Recommended Dilutions	Western Blot 1:200 - 1:1000, Immunohistochemistry 1:50 - 1:200, Immunohistochemistry-Frozen 1:50 - 1:200, Immunocytochemistry 1:50 - 1:200
Application Notes	WB: Detects an approx. 16 kDa protein representing leptin from 3T3-L1 cell extract. IHC: Staining of leptin in sheep brain with NB300-611 results in staining of the ventromedial hypothalamus.




Images

Western Blot: Leptin/OB Antibody [NB300-611] - Analysis of leptin in mouse 3T3-L1.



Immunohistochemistry-Frozen: Leptin/OB Antibody [NB300-611] - Staining of Leptin in sheep brain.



Publications

Kuramoto K, Kim YJ, Hong JH, He C The autophagy protein Becn1 improves insulin sensitivity by promoting adiponectin secretion via exocyst binding Cell reports 2021-05-25 [PMID: 34038729]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-611

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-611

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

