

Product Datasheet

LXR alpha/NR1H3 Antibody

NB300-612

Unit Size: 200 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 10

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-612

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-612



NB300-612**LXR alpha/NR1H3 Antibody**


Product Information	
Unit Size	200 ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS with 1 mg/ml BSA

Product Description	
Description	Novus Biologicals Rabbit LXR alpha/NR1H3 Antibody (NB300-612) is a polyclonal antibody validated for use in IHC, WB, Flow and ICC/IF. Anti-LXR alpha/NR1H3 Antibody: Cited in 7 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	10062
Gene Symbol	NR1H3
Species	Human, Mouse
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 17393442).
Specificity/Sensitivity	Detects recombinant human LXR alpha. This does not detect recombinant human LXR beta.
Immunogen	Synthetic peptide corresponding to residues C D(318) F S Y N R E D F A K A (329) of human LXR alpha.


Product Application Details	
Applications	Western Blot, Flow Cytometry, Gel Super Shift Assays, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:300, Flow Cytometry 3-5ug/10 ⁶ cells, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:50-1:500, Gel Super Shift Assays 1:1 - 1:100
Application Notes	WB: Detects an approx. 64 kDa protein representing recombinant human LXR alpha. IHC usage was reported in the scientific literature (PMID: 17393442). Gel Shift Assays usage was reported in the scientific literature (PMID: 12470667).

Images


Western Blot: LXR alpha/NR1H3 Antibody [NB300-612] - LXRs are expressed in peripheral blood cells. LXR alpha/NR1H3 protein levels in protein extracts from PBMCs from these same donors were detected by Western blotting using rabbit anti-human LXR alpha/NR1H3 polyclonal antisera. Image collected and cropped by CiteAb from the following publication (<https://translational-medicine.biomedcentral.com/articles/10.1186/1479-5876-6-59>), licensed under a CC-BY license.




Immunocytochemistry/Immunofluorescence: LXR alpha/NR1H3 Antibody [NB300-612] - Analysis of LXR alpha (green) showing staining in the nucleus of HeLa cells (right) compared to a negative control without primary antibody (left).




Flow Cytometry: LXR alpha/NR1H3 Antibody [NB300-612] - analysis of LXR alpha was done on HeLa. Cells were fixed with 70% ethanol for 10 minutes, permeabilized with 0.25% Triton X-100 for 20 minutes, and blocked with 2.5% BSA for 30 minutes at room temperature. Cells were labeled with LXR alpha Rabbit Polyclonal Antibody (PA1330, red histogram) or with rabbit isotype control (pink histogram) at 3-5 ug/million cells in 2.5% BSA. After incubation at room temperature for 2 hours, the cells were labeled with Alexa Fluor 488 Goat Anti-Rabbit Secondary Antibody (A11008) at a dilution of 1:400 for 30 minutes at room temperature. The representative 10,000 cells were acquired and analyzed for each sample using an Acoustic Focusing Cytometer. The purple histogram represents unstained control cells and the green histogram represents no-primary-antibody control.




Western Blot: LXR alpha/NR1H3 Antibody [NB300-612] - analysis was performed on nuclear enriched extracts (30 ug lysate) of MCF7 (Lane 1), THP-1 (Lane 2), A549 (Lane 3), SH-SY5Y (Lane 4), HL-60 (Lane 5) and KARPAS 299 (Lane 6). The blot was probed with Rabbit Anti-LXR alpha Polyclonal Antibody (2ug/ml) and detected by chemiluminescence using Goat anti-Rabbit IgG (H+L) Secondary Antibody, HRP conjugate (0.25 ug/ml, 1:4000 dilution). A 50 kDa band corresponding to LXR alpha was observed across the cell lines tested. Known quantity of protein samples was electrophoresed using 4-12 % Bis-Tris gel, Electrophoresis System and Sharp Pre-Stained Protein Standard. Resolved proteins were then transferred onto a nitrocellulose membrane with Blot 2 Dry Blotting System. The membrane was probed with the relevant primary and secondary Antibody following blocking with 5 % skimmed milk. Chemiluminescent detection was performed using ECL Chemiluminescent Substrate Reagent Kit.



Immunocytochemistry/Immunofluorescence: LXR alpha/NR1H3 Antibody [NB300-612] - Analysis of LXR alpha (green) showing staining in the nucleus of MCF-7 cells (right) compared to a negative control without primary antibody (left).



LXRs are expressed in peripheral blood cells. (A) RNA from peripheral blood mononuclear cells obtained from normal human donors was assayed for LXR α and LXR β transcript levels using qPCR. Expression values were normalized to GAPDH levels, represented as the mean \pm SEM. (B) LXR α protein levels in protein extracts from PBMCs from these same donors were detected by Western blotting using rabbit anti-human LXR α polyclonal antisera.



Publications

Rudalska R, Harbig J, Snaebjornsson M et al. LXR alpha activation and Raf inhibition trigger lethal lipotoxicity in liver cancer *Nature Cancer* 2021-02-01 [PMID: 35122079]

Palma GBH, Kaur M miRNA-128 and miRNA-223 regulate cholesterol-mediated drug resistance in breast cancer *IUBMB life* 2023-04-18 [PMID: 37070323] (WB, Human)

Abuirqeba SA The Role of Liver X Receptor in Prostate Cancer Metabolism Thesis 1905-07-13

Bruschi FV, Claudel T, Tardelli M et al. PNPLA3 I148M Variant Impairs Liver X Receptor Signaling and Cholesterol Homeostasis in Human Hepatic Stellate Cells *Hepatol Commun* 2019-09-01 [PMID: 31497741] (WB, Human)

Shen C, Jiang J, Huang C, Zhu W. Gypenoside LVI attenuates foam cell formation by promoting cholesterol export and inhibiting inflammation response *Oncotarget* 2018-10-05

(WB, Mouse)

Poomthavorn P, Wong SHX, Higgins S et al. Activation of a prometastatic gene expression program in hypoxic neuroblastoma cells. *Endocr Relat Cancer*;16(3):991-1004. 2009-01-01 [PMID: 19423615]

Dang H, Liu Y, Pang W et al. Suppression of 2,3-Oxidosqualene Cyclase by High Fat Diet Contributes to Liver X Receptor- α -mediated Improvement of Hepatic Lipid Profile. *J Biol Chem*;284(10):6218-6226. 2009-01-01 [PMID: 19119143]

DiBlasio-Smith EA, Arai M, Quinet EM et al. Discovery and implementation of transcriptional biomarkers of synthetic LXR agonists in peripheral blood cells (BioMed Central). *J Transl Med*. 2008-01-01 [PMID: 18925943]

Zhang, C et al. NO-1886 suppressing atherosclerosis in high-fat-high-sucrose/high-cholesterol fed Bama minipigs is related to upregulating ATP-binding cassette transporter A1. *J Lipid Res*. 2006-01-01 [PMID: 16807312]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-612

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB400-157PEP	LXR alpha/NR1H3 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-612

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

