

Product Datasheet

Phosphodiesterase 4A/PDE4A [p Ser686, p Ser688] Antibody - BSA Free NB300-635-0.05mg

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-635

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-635




NB300-635-0.05mg

Phosphodiesterase 4A/PDE4A [p Ser686, p Ser688] Antibody - BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris/Glycine buffer pH 7.5-7.9, stabilizing protein, and glycerol
Product Description	
Description	Novus Biologicals Rabbit Phosphodiesterase 4A/PDE4A [p Ser686, p Ser688] Antibody - BSA Free (NB300-635) is a polyclonal antibody validated for use in WB, ELISA and IP. Anti-Phosphodiesterase 4A/PDE4A Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5141
Gene Symbol	PDE4A
Species	Human, Mouse, Rat, Rabbit
Specificity/Sensitivity	This detects all known long-form variants of PDE4A, including PDE4A4, PDE4A5, PDE4A8, PDE4Ax when phosphorylated at the PKA phosphorylation site. It will not cross-react with PDE4A phosphorylated at the ERK phosphorylation site. It does not cross-react with PDE4B, PDE4D or any other member of the PDE family.
Immunogen	Synthetic peptide corresponding to PKA phosphorylation site on PDE4A long form enzymes. The peptide was synthesized with phosphorylated serine.
Product Application Details	
Applications	Western Blot, ELISA, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500, ELISA 1:10000, Immunoprecipitation 1:150
Application Notes	Immunoprecipitation and Western Blot: Detects bands at 66 kDa (PDE4A1), 109 kDa (PDE4A5), 106 kDa(PDE4A8), 102 kDa (PDE4Ax) and 76 kDa (testis specific PDE4A variant). Optimal dilutions and concentrations should be determined by the end user.

Images

Western Blot: Phosphodiesterase 4A/PDE4A [p Ser686, p Ser688]
Antibody [NB300-635] - Lane 1: PDE4A Western blot positive Control (PDE4A5); Lane 2: Western blot with Phospho-PDE4A5



Publications

Ale-Agha N, Goy C, Jakobs P et al. CDKN1B/p27 is localized in mitochondria and improves respiration-dependent processes in the cardiovascular system-New mode of action for caffeine PLoS Biol. 2018-06-01 [PMID: 29927970] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-635-0.05mg

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-635

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

