

Product Datasheet

Presenilin-1 Antibody (APS 11) NB300-749

Unit Size: 200ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-749

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-749



NB300-749

Presenilin-1 Antibody (APS 11)


Product Information	
Unit Size	200ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	APS 11
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS with 1 mg/ml BSA

Product Description	
Description	Novus Biologicals Mouse Presenilin-1 Antibody (APS 11) (NB300-749) is a monoclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. Anti-Presenilin-1 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5663
Gene Symbol	PSEN1
Species	Human, Mouse, Rat, Primate
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information.
Specificity/Sensitivity	No cross-reactivity is seen with presenilin 2.
Immunogen	Synthetic peptide corresponding to residues C H(21) L S N T V R S Q N D N R E (34) of human PS1.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:100 - 1:1000, ELISA 4.5 ug/ml, Immunohistochemistry 1:10 - 1:500, Immunocytochemistry/ Immunofluorescence 0.75 ug/ml, Immunoprecipitation 1:10 - 1:500, Immunohistochemistry-Paraffin 0.75 ug/ml
Application Notes	IHC-P: Staining of both the human brain plaque core and dystrophic neurites. WB: Detects an approx. 28 kDa protein representing PS1 N-terminus cleavage product.

Images


Western Blot: Presenilin-1 Antibody (APS 11) [NB300-749] - Analysis of 25 ug of SH-SY5Y (lane 1), human brain (lane 2) and mouse brain (lane 3).




Immunocytochemistry/Immunofluorescence: Presenilin-1 Antibody (APS 11) [NB300-749] - Mouse fibroblasts.



Immunohistochemistry-Paraffin: Presenilin-1 Antibody (APS 11) [NB300-749] - Immunohistochemistry was performed on normal biopsies of deparaffinized Human liver tissue. To expose target proteins, heat induced antigen retrieval was performed using 10mM sodium citrate (pH6.0) buffer, microwaved for 8-15 minutes. Following antigen retrieval tissues were blocked in 3% BSA-PBS for 30 minutes at room temperature. Tissues were then probed at a dilution of 1:200 with a mouse monoclonal antibody recognizing Presenilin 1 or without primary antibody (negative control) overnight at 4C in a humidified chamber. Tissues were washed extensively with PBST and endogenous peroxidase activity was quenched with a peroxidase suppressor. Tissues were counterstained with hematoxylin and prepped for mounting.



Immunohistochemistry-Paraffin: Presenilin-1 Antibody (APS 11) [NB300-749] - Immunohistochemistry was performed on normal biopsies of deparaffinized Human brain tissue. To expose target proteins, heat induced antigen retrieval was performed using 10mM sodium citrate (pH6.0) buffer, microwaved for 8-15 minutes. Following antigen retrieval tissues were blocked in 3% BSA-PBS for 30 minutes at room temperature. Tissues were then probed at a dilution of 1:20 with a mouse monoclonal antibody recognizing Presenilin 1 or without primary antibody (negative control) overnight at 4C in a humidified chamber. Tissues were washed extensively with PBST and endogenous peroxidase activity was quenched with a peroxidase suppressor. Tissues were counterstained with hematoxylin and prepped for mounting.



Publications

Scarpi D, Cirelli D et al. Low molecular weight, non-peptidic agonists of TrkA receptor with NGF-mimetic activity. Cell Death Dis 2012-05-07 [PMID: 22764098] (WB, Rat)

Scarpi D, Cirelli D, Matrone C et al. Low molecular weight, non-peptidic agonists of TrkA receptor with NGF-mimetic activity Cell Death Dis 2012-09-06 [PMID: 22951986] (WB, Rat)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-749

NB820-59177	Human Brain Whole Tissue Lysate (Adult Whole Normal)
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-749

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

