

Product Datasheet

VDAC2 Antibody NB300-871

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-871

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-871



NB300-871

VDAC2 Antibody

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA
Target Molecular Weight	31.6 kDa

Product Description	
Description	Novus Biologicals Goat VDAC2 Antibody (NB300-871) is a polyclonal antibody validated for use in IHC, WB and ELISA. Anti-VDAC2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Goat
Gene ID	7417
Gene Symbol	VDAC2
Species	Human, Rat
Immunogen	Peptide with sequence C-GHKVGLALELEA corresponding to C-Terminus according to NP_003366.2.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Peptide ELISA
Recommended Dilutions	Western Blot 0.3 - 1 ug/mL, Immunohistochemistry, Immunohistochemistry-Paraffin 5 - 6 ug/ml, Peptide ELISA Detection limit 1:2000
Application Notes	WB: Approx 35 kDa band observed in Human Cerebellum lysates and in preliminary testing of CaCo-2 cell lysate (calculated MW of 31.6kDa according to Human NP_003366.2)



Images


Western Blot: VDAC2 Antibody [NB300-871] - Staining (1ug/ml) of Human Cerebellum lysate (35ug protein in RIPA buffer). Detected by chemiluminescence.



Immunohistochemistry-Paraffin: VDAC2 Antibody [NB300-871] - Staining (6ug/ml) of paraffin embedded Rat Liver. Heat induced antigen retrieval with citrate buffer pH 6, HRP-Staining.



Immunohistochemistry-Paraffin: VDAC2 Antibody [NB300-871] - Negative Control showing staining of paraffin embedded Rat Liver, with no primary antibody.



Publications

Hoeffer CA, Dey A, Sachan N et al. The Down Syndrome Critical Region Protein RCAN1 Regulates Long-Term Potentiation Memory via Inhibition of Phosphatase Signaling. J Neurosci 2007-11-28 [PMID: 18045910]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-871

NBL1-17709	VDAC2 Overexpression Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-871

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



