

Product Datasheet

SR-BI/SR-BII Antibody - BSA Free NB400-134

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 11

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB400-134

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB400-134



NB400-134


SR-BI/SR-BII Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Unpurified
Buffer	Whole antisera
Target Molecular Weight	82 kDa
Product Description	
Description	Novus Biologicals Knockout (KO) Validated Rabbit SR-BI/SR-BII Antibody - BSA Free (NB400-134) is a polyclonal antibody validated for use in IHC, WB, Flow, ICC/IF and IP. Anti-SR-BI/SR-BII Antibody: Cited in 11 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	949
Gene Symbol	SCARB1
Species	Human, Mouse, Rat
Reactivity Notes	Rat (PMID: 22097902) and Human (PMID: 17105723) reactivity reported in scientific literature
Immunogen	A peptide from the extracellular domain (residues 230-380) of Scavenger Receptor-BI/BII that was expressed as two tandem copies in bacteria using the pET system.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Block/Neutralize, Immunohistochemistry-Frozen (Negative), Knockout Validated
Recommended Dilutions	Western Blot 1:1000, Flow Cytometry 1:400. Use reported in scientific literature (PMID 22097902), Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:50-1:1000, Immunoprecipitation, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen (Negative) 5 ug/ml, Knockout Validated, Block/Neutralize reported in scientific literature (PMID 24859737)
Application Notes	In Western blot a band is observed at ~ 82 kDa in tissues that express SR-BI and/or SR-BII such as liver, ovary, adrenal glands, and to as lesser extent testes, heart and mammary glands.




Images


Western Blot: SR-BI/SR-BII Antibody [NB400-134] - Detection of RED-1 in 80 ug of total mouse liver lysates. Lane 1: wild-type mice, Lane 2: SR-BI deficient mice.




Immunocytochemistry/Immunofluorescence: SR-BI/SR-BII Antibody [NB400-134] - SR-BI/SR-BII antibody was tested in HeLa cells with DyLight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and Dylight 550 (red).



Immunohistochemistry-Paraffin: SR-BI/SR-BII Antibody [NB400-134] - IHC analysis of formalin-fixed paraffin-embedded tissue section of normal human liver using 5 ug/ml concentration of SR-BI/SR-BII antibody. Specific and expected granular membrane-cytoplasmic staining was observed in the hepatocytes. [40X Magnification]



Flow Cytometry: SR-BI/SR-BII Antibody [NB400-134] - SR-BI/SR-BII antibody was tested at 1:400 in NIH-3T3 cells using an Alexa Fluor 488 secondary (shown in purple). M1 is defined by unstained cells.



Publications

Xiao Wang, Huimin Wang, Bolin Xu, Dong Huang, Chao Nie, Longjun Pu, Gregory J M Zajac, Han Yan, Jingru Zhao, Fangyuan Shi, Brian T Emmer, Jia Lu, Rui Wang, Xiaohui Dong, Jianye Dai, Wenjing Zhou, Chu Wang, Ge Gao, Yan Wang, Cristen Willer, Xiangfeng Lu, Yuangang Zhu, Xiao-Wei Chen Receptor-Mediated ER Export of Lipoproteins Controls Lipid Homeostasis in Mice and Humans. *Cell metabolism* 2021-11-29 [PMID: 33186557]

Kakava S, Schlumpf E, Panteloglou G et al. Brain Endothelial Cells in Contrary to the Aortic Do Not Transport but Degrade Low-Density Lipoproteins via Both LDLR and ALK1 Cells 2022-09-28 [PMID: 36231005] (WB, Human)

May SC, Dron JS, Hegele RA, Sahoo D Human Variant of Scavenger Receptor BI (R174C) Exhibits Impaired Cholesterol Transport Functions *Journal of lipid research* 2021-02-09 [PMID: 33577783]

Emert B, Hasin-Brumshtein Y, Springstead JR et al. HDL inhibits the effects of oxidized phospholipids on endothelial cell gene expression via multiple mechanisms [S] *J Lipid Res.* 2014-01-29 [PMID: 24859737] (BN, Human)

Huang, L;Chambliss, KL;Gao, X;Yuhanna, IS;Behling-Kelly, E;Bergaya, S;Ahmed, M;Michaely, P;Luby-Phelps, K;Darehshouri, A;Xu, L;Fisher, EA;Ge, WP;Mineo, C;Shaul, PW; SR-B1 drives endothelial cell LDL transcytosis via DOCK4 to promote atherosclerosis *Nature* 2019-04-24 [PMID: 31019307] (WB, Human)

Schafer G, Guler R, Murray G et al. The role of scavenger receptor B1 in infection with *Mycobacterium tuberculosis* in a murine model. *PLoS One* 2009-01-01 [PMID: 20041149]

Kartz GA, Holme RL, Nicholson K, Sahoo, D. SR-BI/CD36 Chimeric Receptors Define Extracellular Subdomains of SR-BI Critical for Cholesterol Transport. *Biochemistry* 2014-09-23 [PMID: 25211142]

Norgaard N, Holien T, Jonsson S et al. CpG-oligodeoxynucleotide inhibits Smad-dependent bone morphogenetic protein signaling: effects on myeloma cell apoptosis and in vitro osteoblastogenesis. *J Immunol.* 2010-01-01 [PMID: 20702733]

Hu J, Zhang Z, Shen W-J et al. Differential Roles of Cysteine Residues in Cellular Trafficking, Dimerization, and Function of the HDL Receptor, SR-BI. *Biochemistry.* 2011-11-18 [PMID: 22097902] (FLOW, Rat)

Harder CJ, Meng A, Rippstein P et al. SR-BI undergoes cholesterol-stimulated transcytosis to the bile canaliculus in polarized WIF-B cells. *J Biol Chem* 2007-01-12 [PMID: 17105723] (ICC/IF, Human)

Silver, DL. A carboxyl-terminal PDZ-interacting domain of scavenger receptor B, type I is essential for cell surface expression in liver. *J Biol Chem*;277(37):34042-7. 2002-09-13 [PMID: 12119305] (WB, IP, Human)



Procedures

Western Blot protocol for SR-BI/SR-BII Antibody (NB400-134)

Western Blot Procedure

1. Run 80 ug of protein on a 4-20% Tris-glycine mini-gel at 125V for 60 minutes.
2. Equilibrate gel, nitrocellulose membrane, Whatman paper, and blotting pads in transfer buffer for 15 minutes.
3. Transfer protein to the membrane at 25V for 90 minutes.
4. Allow membrane to air-dry.
5. Block membrane with 1XPBS/5% non-fat milk/0.1% Tween-20 for 1 hour at room temperature (23-27 degrees C).
6. Wash membrane twice, for 5 minutes each, with 1XPBS/0.05% Tween-20 (PBST).
7. Incubate membrane with 1:1,000 dilution of NB400-134 (anti-SR-BI/BII), diluted in 1XPBS/1% BSA, for 1 hour at room temperature.
8. Wash membrane once for 15 minutes, then four times for 5 minutes each, with PBST.
9. Incubate membrane with dilution of goat anti-rabbit IgG-HRP (secondary), diluted in 1XPBS/1% BSA, for 1 hour at room temperature.
10. Wash membrane once for 15 minutes, then four times for 5 minutes each, with PBST.
11. Detect cross-reacting proteins using Chemiluminescence reagents.

Immunocytochemistry/Immunofluorescence protocol for SR-BI/SR-BII Antibody (NB400-134)

Immunocytochemistry Protocol

Culture cells to appropriate density on suitable glass coverslips in 35 mm culture dishes or 6-well plates.

1. Remove culture medium and add 10% formalin to the dish. Fix at room temperature for 5-10 minutes.
2. Remove the formalin and add 0.5% Triton-X 100 in TBS to permeabilize the cells. Incubate for 5-10 minutes.
3. Remove the permeabilization buffer and add wash buffer (i.e. PBS or PBS with 0.1% Tween-20). Be sure to not let the specimen dry out. Gently wash three times for 10 minutes.
4. Alternatively, cells can be fixed with -20C methanol for 10 min at room temperature. Remove the methanol and rehydrate in PBS for 10 min before proceeding.
5. To block nonspecific antibody binding incubate in 10% normal goat serum for 1 hour at room temperature.
6. Add primary antibody at appropriate dilution and incubate at room temperature for 1 hour or at 4 degrees C overnight.
7. Remove primary antibody and replace with wash buffer. Gently wash three times for 10 minutes.
8. Add secondary antibody at the appropriate dilution. Incubate for 1 hour at room temperature.
9. Remove antibody and replace with wash buffer. Gently wash three times for 10 minutes.
10. Nuclei can be staining with 4',6' diamino phenylindole (DAPI) at 0.1 ug/ml, or coverslips can be directly mounted in media containing DAPI.
11. Cells can now be viewed with a fluorescence microscope.

*The above information is only intended as a guide. The researcher should determine what protocol best meets their needs. Please follow proper laboratory procedures for the disposal of formalin.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB400-134

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB400-134

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

