

Product Datasheet

Lck Antibody (LCK-01) - BSA Free NB500-336

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB500-336

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB500-336



NB500-336

Lck Antibody (LCK-01) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	LCK-01
Preservative	15mM Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	Phosphate buffered saline (PBS), pH 7.4
Target Molecular Weight	56 kDa


Product Description	
Description	Novus Biologicals Mouse Lck Antibody (LCK-01) - BSA Free (NB500-336) is a monoclonal antibody validated for use in WB, Flow, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	3932
Gene Symbol	LCK
Species	Human, Mouse (Negative)
Specificity/Sensitivity	The antibody LCK-01 recognizes defined epitope (aa 22-36) of Lck, a 56 kDa Src-family protein tyrosine kinase.
Immunogen	Peptide corresponding to amino acids 22-36 in the sequence of human Lck.

Product Application Details	
Applications	Western Blot, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, CyTOF-ready
Recommended Dilutions	Western Blot 1:1000, Immunocytochemistry/ Immunofluorescence 1:10-1:2000, Immunoprecipitation 1:50, Flow (Intracellular) 1-5 ug/ml, CyTOF-ready
Application Notes	This antibody is Cytof ready.




Images

Western Blot: Lck Antibody (LCK-01) [NB500-336] - Fig. 1. Western Blotting analysis (reducing conditions) of human Lck in whole cell lysate using anti-human Lck (LCK-01). Lane 1: J. CaM-1.6 cell line (a mutant derivate of the JURKAT cell line) transfected with Lck Lane 2: HEK-293T cell line transfected with Lck Lane 3: HEK-293T cell line (non-transfected).



Flow Cytometry: Lck Antibody (LCK-01) [NB500-336] - Staining of Lck in human peripheral blood cells using anti-Lck (LCK-01) Alexa Fluor® 647.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB500-336

NB800-PC2	Jurkat Whole Cell Lysate
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB500-336

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

