

# Product Datasheet

## Glu-Glu Epitope Tag Antibody - BSA Free NB600-354

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

**Publications: 13**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB600-354](http://www.novusbio.com/NB600-354)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB600-354](http://www.novusbio.com/reviews/destination/NB600-354)



**NB600-354****Glu-Glu Epitope Tag Antibody - BSA Free**

<b>Product Information</b>	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Store at 4C. Do not freeze.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.09% Sodium Azide
<b>Isotype</b>	IgG
<b>Purity</b>	Immunogen affinity purified
<b>Buffer</b>	PBS


<b>Product Description</b>	
<b>Description</b>	Novus Biologicals Rabbit Glu-Glu Epitope Tag Antibody - BSA Free (NB600-354) is a polyclonal antibody validated for use in WB, ELISA, ICC/IF and IP. Anti-Glu-Glu Epitope Tag Antibody: Cited in 13 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Rabbit
<b>Species</b>	Epitope Tag
<b>Immunogen</b>	Rabbits were immunized with glu-glu cleavage site (EYMPME) conjugated to KLH. Antibody was isolated by affinity chromatography using the peptide immobilized on solid support. Rabbit anti- human glu-glu affinity purified antibodies were coupled to peptide immobilized on solid support.

<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
<b>Recommended Dilutions</b>	Western Blot 1:1000-1:30000, ELISA Primary-1:1000-1:30000; Coating: 1:100-1:500, Immunocytochemistry/ Immunofluorescence 1:100-1:400, Immunoprecipitation 1 - 4 ug/mg lysate
<b>Application Notes</b>	Suggested working dilutions: * ELISA Coating 1:100-1:500 Primary 1:1000-1:30,000 Western Blots Colorimetric detection 1:1000-1:10,000ECL 1:1000-1:30,000 Cyto 1:100-1:400 *The investigator should determine the optimal working dilution for a specific application. ELISA, Western blot and immunocytochemistry.




## Images


Western Blot: Glu-Glu Epitope Tag Antibody [NB600-354] - Analysis using the HRP conjugate of NB600-354. Detection of 200, 100, or 50 ng of *E. coli* whole cell lysate expressing a multi-tag fusion protein. Antibody used at 0.2 ug/ml (1:5,000).



Western Blot: Glu-Glu Epitope Tag Antibody [NB600-354] - 200, 100, or 50 ng of *E. coli* whole cell lysate expressing a multi-tag fusion protein. Antibody used at 0.04 ug/ml (1:25,000).



Western Blot: Glu-Glu Epitope Tag Antibody [NB600-354] - Triton X-100 (1%0 Whole cell lysate(20mcg) from Cos7 cells transiently transfected with control vector (v) or an expression vector driving expression of Glu-Glu tagged RalBP1 (a to d).



## Publications

Stelmashenko O, Lalo U, Yang Y et al. Activation of trimeric P2X2 receptors by fewer than three ATP molecules. *Mol Pharmacol* 2012-10-01 [PMID: 22828800]

Bradley HJ, Baldwin JM, Goli GR et al. Residues 155 and 348 contribute to the determination of P2X7 receptor function via distinct mechanisms revealed by single-nucleotide polymorphisms. *J Biol Chem* 2011-03-01 [PMID: 21205829]

Roger S, Gillet L, Baroja-Mazo A et al. C-terminal calmodulin-binding motif differentially controls human and rat P2X7 receptor current facilitation. *J Biol Chem* 2010-06-01 [PMID: 20378545]

Li J, Sukumar P, Milligan CJ et al. Interactions, functions, and independence of plasma membrane STIM1 and TRPC1 in vascular smooth muscle cells. *Circ Res* 2008-10-01 [PMID: 18802022]

Xia R, Mei ZZ, Mao HJ et al. Identification of pore residues engaged in determining divalent cationic permeation in transient receptor potential melastatin subtype channel 2. *J Biol Chem* 2008-10-01 [PMID: 18687688]

Roger S, Pelegrin P, Surprenant A et al. Facilitation of P2X7 receptor currents and membrane blebbing via constitutive and dynamic calmodulin binding. *J Neurosci* 2008-06-01 [PMID: 18562610]

Liu X, Surprenant A, Mao HJ et al. Identification of key residues coordinating functional inhibition of P2X7 receptors by zinc and copper. *Mol Pharmacol* 2008-01-01 [PMID: 17959713]

Guo C, Masin M, Qureshi OS et al. Evidence for functional P2X4/P2X7 heteromeric receptors. *Mol Pharmacol* 2007-12-01 [PMID: 17785580]

Mei ZZ, Xia R, Beech DJ et al. Intracellular coiled-coil domain engaged in subunit interaction and assembly of melastatin-related transient receptor potential channel 2. *J Biol Chem* 2006-12-01 [PMID: 17060318]

Young MT, Pelegrin P, Surprenant A et al. Amino acid residues in the P2X7 receptor that mediate differential sensitivity to ATP and BzATP. *Mol Pharmacol* 2007-01-01 [PMID: 17032903]

Young MT, Pelegrin P, Surprenant A et al. Identification of Thr283 as a key determinant of P2X7 receptor function. *Br J Pharmacol* 2006-10-01 [PMID: 16940988]

Mei ZZ, Mao HJ, Jiang LH et al. Conserved cysteine residues in the pore region are obligatory for human TRPM2 channel function. *Am J Physiol Cell Physiol* 2006-11-01 [PMID: 16822940]

More publications at <http://www.novusbio.com/NB600-354>





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB600-354**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB600-354](http://www.novusbio.com/reviews/submit/NB600-354)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

