

Product Datasheet

AU1 Epitope Tag Antibody - BSA Free NB600-453

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 15

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB600-453

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB600-453



NB600-453**AU1 Epitope Tag Antibody - BSA Free**

Product Information	
Unit Size	0.1 ml
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit AU1 Epitope Tag Antibody - BSA Free (NB600-453) is a polyclonal antibody validated for use in WB, ELISA, ICC/IF and IP. Anti-AU1 Epitope Tag Antibody: Cited in 15 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Species	Epitope Tag
Immunogen	Rabbits were immunized with AU1 (DTYRYI) conjugated to KLH. Antibody was isolated by affinity chromatography using the peptide immobilized on solid support.

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000 - 1:3000, ELISA 1:1000 - 1:3000, Immunocytochemistry/ Immunofluorescence 1:100 - 1:400, Immunoprecipitation 1 - 4 ug/mg lysate
Application Notes	AU1 Epitope Tag antibody validated for WB, IF from a verified customer reviews.




Images

Western Blot: AU1 Epitope Tag Antibody [NB600-453] - 200, 100, or 50 ng of *E. coli* whole cell lysate expressing a multi-tag fusion protein. Antibody used at 0.04 ug/mL (1:25,000).



Western Blot: AU1 Epitope Tag Antibody [NB600-453] - analysis of overexpressed AU1-tagged recombinant proteins in S2 cell extract. Expected lanes are marked with asterisk (*). WB image submitted by a verified customer review.



Publications

Petersen M, Lotke R, Hopfensperger K et al. Inhibition of Infectious HIV-1 Production by Rerouting the Cellular Furin Inhibitor Serpin B8 *Journal of virology* 2023-06-29 [PMID: 37272794]

Ikegame S, Carmichael JC, Wells H et al. Metagenomics-enabled reverse-genetics assembly and characterization of myotis bat morbillivirus *Nature microbiology* 2023-06-01 [PMID: 37142773] (WB, Epitope Tag)

Details:

Dilutions: Flow 1:1,000, WB: 1:1,000

Lee B, Ikegame S, Carmichael J Et al. Zoonotic potential of a novel bat morbillivirus *Research square* 2021-09-29 [PMID: 34611656] (WB)

Pryce R, Azarm K et al. A key region of molecular specificity orchestrates unique ephrin-B1 utilization by Cedar virus. *Life Sci Alliance* 2020-01-01 [PMID: 31862858] (WB)

Hotter D, Bosso M, Jonsson KL et al. IFI16 Targets the Transcription Factor Sp1 to Suppress HIV-1 Transcription and Latency Reactivation *Cell Host Microbe* 2019-06-12 [PMID: 31175045]

Zhang T, Wu C, Yang X et al. Guanylate-Binding Proteins 2 and 5 Exert Broad Antiviral Activity by Inhibiting Furin-Mediated Processing of Viral Envelope Proteins *Cell Rep* 2019-05-14 [PMID: 31091448] (WB, Human)

Heusinger E, Deppe k, Sette P et al. Preadaptation of SIVsmm facilitated Env-mediated counteraction of human tetherin by HIV-2. *J. Virol.* 2018-07-05 [PMID: 29976668] (WB, Human)

Hotter D, Krabbe T, Reith E et al. Primate lentiviruses use at least three alternative strategies to suppress NF- κ B-mediated immune activation *PLoS Pathog.* 2017-08-01 [PMID: 28859166] (WB, Human)

Mack K, Starz K, Sauter D et al. Efficient Vpu-mediated tetherin antagonism by an HIV-1 group O strain. *J Virol* 2017-01-11 [PMID: 28077643]

HirmondoR, Szabo JE, Nyiri K et al. Cross-species inhibition of dUTPase via the Staphylococcal Stl protein perturbs dNTP pool and colony formation in *Mycobacterium* DNA Repair (Amst.) 2015-03-19 [PMID: 25841100] (WB)

Rona G, Palinkas HI, Borsos M et al. NLS copy-number variation governs efficiency of nuclear import - case study on dUTPases. *FEBS J.* 2014-10-04 [PMID: 25283549] (ICC/IF)

Kingston D, Chang H, Ensser A et al. Inhibition of retromer activity by herpesvirus saimiri tip leads to CD4 downregulation and efficient T cell transformation. *J Virol* 2011-10-01 [PMID: 21849449]

More publications at <http://www.novusbio.com/NB600-453>



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB600-453

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB600-453

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

