

Product Datasheet

UbcH8/Ube2L6 Antibody (k1H3) - BSA Free NBP1-04351

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-04351

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-04351



NBP1-04351

UbcH8/Ube2L6 Antibody (k1H3) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	k1H3
Preservative	0.02% Sodium Azide
Isotype	IgG2b Kappa
Purity	Protein G purified
Buffer	PBS (pH 7.4), 10% Glycerol

Product Description	
Description	Novus Biologicals Mouse UbcH8/Ube2L6 Antibody (k1H3) - BSA Free (NBP1-04351) is a monoclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9246
Gene Symbol	UBE2L6
Species	Human
Immunogen	Recombinant human UbcH8/Ube2L6 (1-152aa) purified from E. coli


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000, ELISA 1:100-1:2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:10-1:500

Images


Western Blot: UbcH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - The Cell lysate (20ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human UBE2L6 antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: Jurkat cell lysate




Immunocytochemistry/Immunofluorescence: Ubch8/Ube2L6 Antibody (k1H3) [NBP1-04351] - ICC/IF analysis of UBE2L6 in HeLa cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human UBE2L6 antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).




Immunohistochemistry-Paraffin: Ubch8/Ube2L6 Antibody (k1H3) [NBP1-04351] - Ube2L6 Antibody (k1H3) [NBP1-04351]-Paraffin embedded sections of human breast lobule tissue were incubated with anti-human UBE2L6 for 2 hours at room temperature. Antigen retrieval was performed in 0.1M sodium citrate buffer and detected using Diaminobenzidine (DAB)



Western Blot: Ubch8/Ube2L6 Antibody (k1H3) [NBP1-04351] - Ube2L6 Antibody (k1H3) [NBP1-04351]-Cell lysates of MCF7 (20ug) were resolved by SDS-PAGE, transferred to NC membrane and probed with anti-human UBE2L6 . Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.



Western Blot: Ubch8/Ube2L6 Antibody (k1H3) [NBP1-04351] - The Cell lysates (5ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human UBE2L6 antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: 293T cell lysate Lane 2.: UBE2L6 Transfected 293T cell lysate





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-04351

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-04351

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

