

Product Datasheet

Pdi1p Antibody (38H8) NBP1-05391

Unit Size: 0.5 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-05391

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-05391



NBP1-05391

Pdi1p Antibody (38H8)

Product Information	
Unit Size	0.5 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	38H8
Preservative	0.035% Sodium Azide
Isotype	IgG1
Purity	Tissue culture supernatant
Buffer	Supplied as concentrated hybridoma cell culture media.
Target Molecular Weight	58.2 kDa

Product Description	
Description	Novus Biologicals Mouse Pdi1p Antibody (38H8) (NBP1-05391) is a monoclonal antibody validated for use in WB, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Species	Yeast
Reactivity Notes	<i>S. cerevisiae</i>
Marker	Endoplasmic Reticulum Marker
Immunogen	Yeast nuclear preparations

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:10000, Immunocytochemistry/ Immunofluorescence 1:500-1:1000, Immunoprecipitation
Application Notes	This Pdi1p (38H8) antibody is useful in Western blot and Immunocytochemistry/Immunofluorescence. In WB of yeast protein extracts a ~65 kDa band can be seen. In ICC/IF using yeast cells, endoplasmic reticulum staining can be seen.

Images

Western Blot: Pdi1p Antibody (38H8) [NBP1-05391] - Western Blot analysis of whole yeast protein extracts detected Pdi1p. Molecular weight detected is 58.2 kDa.

38H8
Pdi1p
(58.2)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-05391

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-05391

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

