

Product Datasheet

HPI4 Antibody (HIC1-5F10) - BSA Free NBP1-18948

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-18948

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-18948



NBP1-18948


HPi4 Antibody (HIC1-5F10) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	HIC1-5F10
Preservative	0.05% Sodium Azide
Isotype	IgG1
Purity	Tissue culture supernatant
Buffer	Tissue culture supernatant
Product Description	
Description	Novus Biologicals Mouse HPi4 Antibody (HIC1-5F10) - BSA Free (NBP1-18948) is a monoclonal antibody validated for use in IHC, Flow and ICC/IF. Anti-HPi4 Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Species	Human, Mouse
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 26515645).
Marker	Endocrine pan-islet cell Marker
Specificity/Sensitivity	Specific for multiple endocrine cell types.
Immunogen	Human pancreatic enriched islet cells containing low levels of exocrine and ductal cells.
Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Flow Cytometry 1:50-1:100, Immunohistochemistry 1:100, Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunohistochemistry-Frozen 1:100
Application Notes	This HPi4 (HIC1-5F10) antibody is useful for Immunohistochemistry on acetone fixed frozen sections, Immunocytochemistry/Immunofluorescence and Flow cytometry.




Images

Immunocytochemistry/Immunofluorescence: HPi4 Antibody (HIC1-5F10) [NBP1-18948] - Immunofluorescence on frozen section of human pancreas.



Flow Cytometry: HPi4 Antibody (HIC1-5F10) [NBP1-18948] - Analysis of enzyme dispersed human pancreas cells.



Publications

Korytnikov R, Nostro MC. Generation of polyhormonal and multipotent pancreatic progenitor lineages from human pluripotent stem cells. *Methods*. 2015-10-27 [PMID: 26515645]

Dorrell, et al. Isolation of major pancreatic cell types long-term culture-initiating cells using novel human surface markers. *Stem Cell Rsch* 1(3):183-194. 2008-01-01 [PMID: 19383399] (FLOW, IHC-Fr, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-18948

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB7539	Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-18948G	HPi4 Antibody (HIC1-5F10) [DyLight 488]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-18948

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

