

Product Datasheet

WDR5 Antibody - BSA Free

NBP1-30063

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-30063

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-30063



NBP1-30063

WDR5 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)


Product Description	
Description	Novus Biologicals Rabbit WDR5 Antibody - BSA Free (NBP1-30063) is a polyclonal antibody validated for use in IHC and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	11091
Gene Symbol	WDR5
Species	Human
Immunogen	maps to a region between residue 1 and 50 of human WD repeat domain 5 using the numbering given in entry NP_060058.1

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation, Western Blot (Negative)
Recommended Dilutions	Immunohistochemistry 1:500 - 1:2000, Immunoprecipitation 2 - 10 ug/mg lysate, Immunohistochemistry-Paraffin 1:500 - 1:2000, Western Blot (Negative)
Application Notes	Epitope retrieval with Tris-EDTA pH 9.0 is recommended for FFPE tissue sections. Western Blot Not recommended.




Images

Immunohistochemistry-Paraffin: WDR5 Antibody [NBP1-30063] -
 Sample: FFPE section of human ovarian carcinoma. Antibody: Affinity
 purified rabbit anti-fs18 WDR5 fs20 used at a dilution of 1:1,000 (1ug/ml).
 Detection: DAB



Immunoprecipitation: WDR5 Antibody [NBP1-30063] - Detection of
 human WDR5 by western blot of immunoprecipitates. Samples: Whole
 cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from HeLa
 cells. Antibodies: Affinity purified rabbit anti-WDR5 antibody A302-429A
 (lot A302-429A-2) used for IP at 6 ug per reaction. WDR5 was also
 immunoprecipitated by rabbit anti-WDR5 antibody A302-430A. For
 blotting immunoprecipitated WDR5, A302-430A was used at 1 ug/ml.
 Detection: Chemiluminescence with an exposure time of 3 minutes.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-30063

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-30063

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

