

Product Datasheet

Kv7.1 Antibody NBP1-30119

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-30119

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-30119



NBP1-30119**Kv7.1 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA

Product Description	
Description	Novus Biologicals Goat Kv7.1 Antibody (NBP1-30119) is a polyclonal antibody validated for use in IHC, WB and ELISA. Anti-Kv7.1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Goat
Gene ID	3784
Gene Symbol	KCNQ1
Species	Human
Specificity/Sensitivity	This antibody is expected to recognize both reported isoforms (NP_000209.2; NP_861463.1).
Immunogen	Peptide with sequence C-EQLTVPRRGPDGECS corresponding to C-Terminus according to NP_861463.1.

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Peptide ELISA, Control
Recommended Dilutions	Western Blot 1.0 - 3.0 ug/ml, Immunohistochemistry 6.0 - 8.0 ug/ml, Peptide ELISA Detection limit 1:16000, Control
Application Notes	WB: Approx. 60 kDa band observed in human heart lysates (calculated MW of 61.5 kDa band according to NP_861463.1). An additional band of unknown identity was also consistently observed at 28 kDa. This band was successfully blocked by incubation with the immunizing peptide.



Images


Western Blot: Kv7.1 Antibody [NBP1-30119] - Analysis of Kv7.1 in Human Heart lysate (35ug protein in RIPA buffer) using NBP1-30119 at 1 ug/ml. Primary incubation was 1 hour. Detected by chemiluminescence.



Immunohistochemistry: Kv7.1 Antibody [NBP1-30119] - (8ug/ml) staining of paraffin embedded Rat Liver. Heat induced antigen retrieval with citrate buffer Ph 6, HRP-staining.



Control: Kv7.1 Antibody [NBP1-30119] - Negative Control showing staining of paraffin embedded Rat Liver, with no primary antibody.



Publications

Yasuda K, Miyake K, Horikawa Y et al. Variants in KCNQ1 are associated with susceptibility to type 2 diabetes mellitus. Nat Genet 2008-08-17 [PMID: 18711367]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-30119

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-30119

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



