

# Product Datasheet

## Hexokinase 1 Antibody - BSA Free NBP1-31600

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-31600](http://www.novusbio.com/NBP1-31600)

Updated 9/25/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-31600](http://www.novusbio.com/reviews/destination/NBP1-31600)



**NBP1-31600**


Hexokinase 1 Antibody - BSA Free

Product Information	
<b>Unit Size</b>	100 ul
<b>Concentration</b>	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.01% Thimerosal
<b>Isotype</b>	IgG
<b>Purity</b>	Antigen Affinity-purified
<b>Buffer</b>	0.1M Tris, 0.1M Glycine, 10% Glycerol
<b>Target Molecular Weight</b>	102 kDa
Product Description	
<b>Description</b>	Novus Biologicals Rabbit Hexokinase 1 Antibody - BSA Free (NBP1-31600) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Hexokinase 1 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Rabbit
<b>Gene ID</b>	3098
<b>Gene Symbol</b>	HK1
<b>Species</b>	Human, Zebrafish
<b>Reactivity Notes</b>	Bovine reactivity reported in scientific literature (PMID: 31261667). Immunogen displays the following percentage of sequence identity for non-tested species: Chicken (88%), Xenopus laevis (86%).
<b>Immunogen</b>	Recombinant protein encompassing a sequence within the center region of human Hexokinase 1. The exact sequence is proprietary.
Product Application Details	
<b>Applications</b>	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry Whole-Mount
<b>Recommended Dilutions</b>	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Immunohistochemistry Whole-Mount
<b>Application Notes</b>	IHC-WM-Assay dependent




## Images


Western Blot: Hexokinase 1 Antibody [NBP1-31600] - Sample (30 ug of whole cell lysate) C: zebrafish muscle 7.5% SDS PAGE diluted at 1:1000



Immunocytochemistry/Immunofluorescence: Hexokinase 1 Antibody [NBP1-31600] - Analysis of methanol-fixed HeLa, using HKX I antibody (Green) at 1:500 dilution. Alpha-tubulin filaments were labeled with an alpha Tubulin antibody (Red) at 1:2000.



Immunohistochemistry Whole-Mount: Hexokinase 1 Antibody [NBP1-31600] - Hexokinase 1 antibody detects Hexokinase 1 protein on zebrafish by whole mount immunohistochemical analysis. Sample: 2 days-post-fertilization zebrafish embryo. Hexokinase 1 antibody dilution: 1:100.




Immunohistochemistry-Paraffin: Hexokinase 1 Antibody [NBP1-31600] - Human breast cancer, using HKX I antibody at 1:250 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.



Immunohistochemistry-Paraffin: Hexokinase 1 Antibody [NBP1-31600] - Analysis of paraffin-embedded zebrafish tissue, using Hexokinase 1 antibody at 1:300 dilution.



Western Blot: Hexokinase 1 Antibody [NBP1-31600] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with Hexokinase 1 antibody diluted at 1:3000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.



## Publications

Rajamanickam GD, Kastelic JP, Thundathil JC Testis-Specific Isoform of Na/K-ATPase (ATP1A4) Interactome in Raft and Non-Raft Membrane Fractions from Capacitated Bovine Sperm Int J Mol Sci 2019-06-28 [PMID: 31261667] (WB, Bovine)



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-31600**

---

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-31600](http://www.novusbio.com/reviews/submit/NBP1-31600)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

