

Product Datasheet

hnRNP H Antibody - BSA Free NBP1-31648

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31648

Updated 9/25/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31648



NBP1-31648

hnRNP H Antibody - BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Target Molecular Weight	49 kDa

Product Description	
Description	Novus Biologicals Rabbit hnRNP H Antibody - BSA Free (NBP1-31648) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-hnRNP H Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	3187
Gene Symbol	HNRNPH1
Species	Human, Mouse
Reactivity Notes	Xenopus laevis (91%), Zebrafish (87%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human hnRNP H. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Gel Supershift Assay
Recommended Dilutions	Western Blot 1:5000-1:20000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Gel Supershift Assay
Application Notes	Use in Gel super shift assay reported in scientific literature (PMID: 30503283).

Images


Western Blot: hnRNP H Antibody [NBP1-31648] - Various whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with hnRNP H antibody diluted by 1:1000.




Immunocytochemistry/Immunofluorescence: hnRNP H Antibody [NBP1-31648] - HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: hnRNP H protein stained by hnRNP H antibody [N1C1] diluted at 1:200. Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [B-5-1-2] diluted at 1:10000. Blue: Hoechst 33342 staining.




Immunohistochemistry-Paraffin: hnRNP H Antibody [NBP1-31648] - Paraffin-embedded NCIN87 xenograft, using antibody at 1:100 dilution.




Western Blot: hnRNP H Antibody [NBP1-31648] - Sample (50 ug of whole cell lysate) A: Mouse brain 10% SDS PAGE; antibody diluted at 1:1000.




Western Blot: hnRNP H Antibody [NBP1-31648] - Sample (30 ug of whole cell lysate) A: A431 B: H1299 C: Hela D: Hep G2 10% SDS PAGE; antibody diluted at 1:1000.



Immunocytochemistry/Immunofluorescence: hnRNP H Antibody [NBP1-31648] - A431 cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: hnRNP H protein stained by hnRNP H antibody [N1C1] diluted at 1:500. Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [114] (NBP2-43837) diluted at 1:1000. Blue: Hoechst 33342 staining.



Western Blot: hnRNP H Antibody [NBP1-31648] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with hnRNP H antibody [N1C1] diluted at 1:10000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Publications

Wang ZF, Ursu A, Childs-Disney JL et al. The Hairpin Form of r(G4C2)_{exp} in c9ALS/FTD Is Repeat-Associated Non-ATG Translated and a Target for Bioactive Small Molecules Cell Chem Biol 2018-11-29 [PMID: 30503283] (GS, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-31648

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31648

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

