

Product Datasheet

PDK3 Antibody NBP1-32581

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 5

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32581

Updated 9/25/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32581



NBP1-32581**PDK3 Antibody**


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol
Target Molecular Weight	47 kDa

Product Description	
Description	Novus Biologicals Rabbit PDK3 Antibody (NBP1-32581) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-PDK3 Antibody: Cited in 5 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5165
Gene Symbol	PDK3
Species	Human, Mouse, Squirrel, Xenopus
Reactivity Notes	Frog (100%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human PDK3. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: PDK3 Antibody [NBP1-32581] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with PDK3 antibody [N1C2] diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunohistochemistry-Paraffin: PDK3 Antibody [NBP1-32581] - Ca922 xenograft, using PDK3 antibody at 1:100 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



Western Blot: PDK3 Antibody [NBP1-32581] - Mouse tissue extract (50 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with PDK3 antibody [N1C2] diluted at 1:500.



Western Blot: PDK3 Antibody [NBP1-32581] - Whole cell extract (30 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with PDK3 antibody [N1C2] diluted at 1:1000.



Immunofluorescence analysis of methanol-fixed HepG2, using PDK3 (NBP1-32581) antibody at 1:200 dilution.



Publications

□korja Mili? N, Dolinar K, Mi□ K et al. Suppression of Pyruvate Dehydrogenase Kinase by Dichloroacetate in Cancer and Skeletal Muscle Cells Is Isoform Specific and Partially Independent of HIF-1? International Journal of Molecular Sciences 2021-08-10 [PMID: 34445316] (Human)

Yang ES, Do Y, Cheon SY et al. Andrographolide suppresses aerobic glycolysis and induces apoptotic cell death by inhibiting pyruvate dehydrogenase kinase 1 expression Oncology reports 2023-04-01 [PMID: 36825595] (WB, Human)

Details:

Dilution used in WB 1:000 in 5% skim milk

Schoonjans CA, Joudiou N, Brusa D et al. Acidosis-induced metabolic reprogramming in tumor cells enhances the anti-proliferative activity of the PDK inhibitor dichloroacetate Cancer Lett. 2020-02-01 [PMID: 31812695] (WB, Human)

Hou X, Zhang L, Han L et al. Differential function of pyruvate dehydrogenase kinases during mouse oocyte maturation. J. Cell. Sci. 2015-05-19 [PMID: 25991547] (ICC/IF, Mouse)

Ho N, Coomber BL. Pyruvate dehydrogenase kinase expression and metabolic changes following dichloroacetate exposure in anoxic human colorectal cancer cells. Exp Cell Res 2015-02-01 [PMID: 25536473] (Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32581

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32581

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

