

Product Datasheet

P2X7/P2RX7 Antibody NBP1-37775

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-37775

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-37775



NBP1-37775**P2X7/P2RX7 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA


Product Description	
Description	Novus Biologicals Goat P2X7/P2RX7 Antibody (NBP1-37775) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. Anti-P2X7/P2RX7 Antibody: Cited in 4 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Goat
Gene ID	5027
Gene Symbol	P2RX7
Species	Human, Mouse
Reactivity Notes	Use in Mouse reported in scientific literature (PMID:34391810).
Immunogen	Peptide with sequence YETNKVTRIQSMNY-C corresponding to N-Terminus according to NP_002553.2.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Peptide ELISA
Recommended Dilutions	Western Blot 1 - 3 ug/mL, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin 3.75 ug/mL, Peptide ELISA detection limit 1:16000
Application Notes	Use in ICC/IF reported in scientific literature (PMID:34391810)




Images


Western Blot: P2X7/P2RX7 Antibody [NBP1-37775] - Staining of Human Brain (Frontal Cortex) lysate (35 ug protein in RIPA buffer) with (B) and without (A) blocking with the immunizing peptide. Antibody at 0.3 ug/mL. Primary incubation was 1 hour. Detected by chemiluminescence.



Immunohistochemistry-Paraffin: P2X7/P2RX7 Antibody [NBP1-37775] - Staining of paraffin embedded Human Kidney. Antibody at 3.75 ug/mL. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.



Immunohistochemistry-Paraffin: P2X7/P2RX7 Antibody [NBP1-37775] - Staining of paraffin embedded Human Brain Cortex. Antibody at 3.75 ug/mL. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.



Publications

Alnafisah R, Lundh A, Asah SM et al. Altered purinergic receptor expression in the frontal cortex in schizophrenia Schizophrenia (Heidelberg, Germany) 2022-11-14 [PMID: 36376358] (WB, Human)

Details:

Dilution used in WB 1:1000

Carvalho K, Martin E, Ces A Et al. P2X7-deficiency improves plasticity and cognitive abilities in a mouse model of Tauopathy Progress in neurobiology 2021-08-12 [PMID: 34391810] (ICC/IF, IF/IHC, Mouse)

Martin E, Amar M, Dalle C, et al New role of P2X7 receptor in an Alzheimer's disease mouse model. Mol Psychiatry 2018-06-22 [PMID: 29934546] (IHC-P, Human)

Tomasinsig L, Pizzirani C, Skerlavaj B et al. The human cathelicidin LL-37 modulates the activities of the P2X7 receptor in a structure-dependent manner. J Biol Chem 2008-11-01 [PMID: 18765670]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-37775

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-37775

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



